

Protocol for use No. 4: upcyte® Hepatocyte Cytotoxicity Testing

Introduction

This PFU describes how to perform a cytotoxicity assay using upcyte® Hepatocytes cultured in 96-well plates.

Required products

- **upcyte Hepatocytes (cryopreserved)**

Each vial contains $\sim 5 \times 10^6$ upcyte Hepatocytes (upcyte; "upregulated hepatocytes") frozen per vial from which at least 50% recovery is expected after thawing.

- **Hepatocyte Thawing Medium**

This is ready-to-use for thawing upcyte Hepatocytes. No additional supplements are required.

- **Hepatocyte High Performance Medium**

The Hepatocyte High Performance Medium is designed for the optimal culture and incubation of compounds with upcyte Hepatocytes. In order to obtain Hepatocyte High Performance Medium add the entire contents of supplement A and L-Glutamine to the basal medium.

Due to the manufacturing process the medium may appear opaque but this does not affect the performance of the cells.

The shelf life of the fully supplemented medium is 6 weeks. Once fully supplemented do not freeze the Hepatocyte High Performance Medium. Add antibiotics if preferred.

- **Collagen coated (Type I) culture plates**

These are not provided by BioIVT. Either use pre-coated plates or prepare them as follows: dilute collagen type I with 0.02M acetic acid to a final concentration of 50µg/mL. Add 0.1mL/cm² of the diluted collagen solution to the culture dishes and incubate for 1h at room temperature. Wash the plate twice with PBS and use directly or air dry before storing at 4°C.

Storage conditions

upcyte Hepatocytes can be stored in liquid or vapour phase nitrogen. They should not be stored at -70°C.

Store basal medium as well as fully supplemented Hepatocyte High Performance Medium protected from light at 2 – 8°C. Store Supplement A at -20°C. The expiration date is indicated on the label of the basal medium as well as on the supplement label.

Thawing of cryopreserved upcyte Hepatocytes

1. Thaw and seed upcyte hepatocytes according to [PFU No. 12](#).
2. Incubate cells for 7 days with medium change according to [PFU No. 12](#).

Treatment of upcyte Hepatocytes for cytotoxicity screening

1. Remove the medium and replace with 0.15mL Hepatocyte High Performance Medium containing either vehicle control or test compound.
2. We recommend incubating treated cells for 72-96h.
3. Determine the viability of cells using standard endpoint measurements (e.g. ATP, LDH, MTT/MTS).

Results using upcyte Hepatocytes





Figure 1 shows an example of the cytotoxicity of chlorpromazine, furazolidone and aspirin in upcyte Hepatocytes from Donor 422A. This shows that these cells can discriminate between toxic and non-toxic compounds.



The cytotoxicity of chlorpromazine in different upcyte Hepatocyte donors is shown in Figure 2. The EC₅₀ for this compound using ATP as the endpoint is between 16 – 30 µM.

Product information

Product	Supplements/Components	Product number
upcyte Hepatocytes	<ul style="list-style-type: none"> cryopreserved (1 mL) 	CHE002
Hepatocyte Thawing Medium	<ul style="list-style-type: none"> ready-to-use (50mL) 	MHE001
Hepatocyte High Performance Medium	<ul style="list-style-type: none"> Basal Medium (500mL) Supplement A (5mL) L-Glutamine (5 mL) 	MHE003

Purchaser Notification

Limited Use Label License (upcyte technologies)

This cellular product (upcyte® cells) generated using the upcyte technology is provided under an intellectual property license from upcyte technologies GmbH. The transfer of the upcyte cells is conditioned on the buyer using the purchased product solely in research conducted by the buyer, and the buyer must not use, sell or otherwise transfer this product, its components or any data generated with the product

(a) for diagnostic, therapeutic or prophylactic purposes; or

(b) for resale; or

(c) for the production of therapeutic, diagnostic, prophylactic or any other products; or

(d) to provide a service to deliver information or materials to a third party.

This intellectual property license is timely limited and granted to the buyer for 12 month after delivery of the upcyte cells. After this time the buyer has to destroy remaining upcyte cells and populations thereof or obtain a follow-up license. For information on purchasing a license to this upcyte technologies product for other purposes, for an extended time frame and/or for commercial purposes, contact upcyte technologies GmbH, Osterfeldstrasse 12-14, 22529 Hamburg, Germany or info@upcyte.com.

During the generation of these upcyte cells a product licensed from Life Technologies Corporation containing the following Limited Use Label License was used:

Limited Use Label License (Life Technologies Corporation)

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, and the buyer must not use, sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; or (b) resale, whether or not resold for use in research or (c) or for the commercial production of therapeutic, diagnostic, or prophylactic products or other products not intended for research use; or (d) to provide a service to deliver information or materials to a customer other than for research purposes of such customer. For information on purchasing a license to this product for purposes other than research, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.

