Use and Care Instructions

for your Invisalign® aligners

ENGLISH.....	4	ITALIANO.....	77
العربية.....	9	LATVISKI.....	82
БЪЛГАРСКИ.....	13	LIETUVIŲ K.	86
HRVATSKI.....	18	NORSK.....	91
CATALÀ.....	23	POLSKI.....	96
ČESKY.....	28	PORTUGUÊS.....	101
DANSK.....	32	ROMÂNĂ.....	106
NEDERLANDS.....	37	РУССКИЙ.....	110
EESTI.....	42	SRPSKI.....	116
SUOMI.....	47	SLOVENČINA.....	120
FRANÇAIS.....	51	SLOVENŠČINA.....	125
DEUTSCH.....	57	ESPAÑOL.....	130
ΕΛΛΗΝΙΚΑ.....	62	SVENSKA.....	135
עברית.....	69	TÜRKÇE.....	140
MAGYAR.....	72	УКРАЇНСЬКА.....	144


Figure 1


Figure 2


Figure 3


Figure 8


Figure 9


Figure 10


Figure 11


Figure 4


Figure 5


Figure 6


Figure 7

Figure 12

Figure 13

Figure 14

ENGLISH

Read these instructions before use

Before using Invisalign® aligners, patients (and, when applicable, their parents or legal guardians) should carefully read these Use and care instructions for important information about aligner wear including warnings and precautions, proper care and maintenance, and oral hygiene requirements during treatment. We recommend that you save these instructions for future reference. If you have questions or concerns, please contact your doctor.

What is the Invisalign® system?

The Invisalign® System consists of a series of doctor-prescribed, thin, clear plastic removable orthodontic appliances (aligners) that gently move the patient's teeth in small increments from their original state to a more optimal, treated state to address malocclusion. The first device in the series matches the patient's dentition in its current state and then each subsequent aligner stage has the shape of the dentition shifted gradually toward the final desired position. Each aligner covers the teeth and the gingival tissue between the teeth. The aligners are trimmed near the gum line between the teeth and the gingival tissue so that the edge of each appliance will be minimally visible. Each plastic aligner surrounds the teeth captured by the prescribing dentist in their impression/ intraoral scan, and the aligners may also touch the gingival tissue near the teeth of the patient. The aligners can accommodate use of tooth attachments and elastics (through use of precision cut outs or buttons or power arms). The system can also be ordered with additional aligner features such as precision cut outs, bite ramps, compliance indicators and mandibular advancement features, such as precision wings and occlusal blocks, position the patient's jaw to address skeletal malocclusion. The choice of precision wings or occlusal blocks for mandibular advancement is made based on the doctor's discretion. The system with Invisalign aligners is used in patients with primary, mixed (primary and permanent), and permanent dentition. Invisalign aligners with mandibular advancement feature(s) are used in patients with mixed and permanent dentition.

For patients who are minors, your parent or guardian or another responsible adult should review these instructions and help you follow these directions.

Intended Purpose/Intended User/ Patient population/Intended use environment:

Intended Purpose:

The Invisalign System is intended for the orthodontic treatment of malocclusion. The optional mandibular advancement feature(s) are indicated to treat skeletal malocclusion.

Intended Patient population:

The system is intended to be used in patients with primary, mixed (primary and permanent), and permanent dentition. Invisalign aligners with mandibular advancement feature(s) are used in patients with mixed and permanent dentition. The Invisalign system should not be used in patients with active periodontal disease.

Intended Users:

Children, Adolescents & Adults with malocclusion.

Intended use environment:

After initial distribution in the dental professional's office, the Invisalign System is intended for home use.

Performance characteristics:

Orthodontic movement of teeth of up to 0.25 mm per aligner.

Orthodontic treatment of malocclusion

Clinical Benefit:

An improvement in dental alignment is the clinical attribute that leads to the more subjective benefits of orthodontic treatment promoted by the orthodontic specialty such as a beautiful smile, improved self-esteem and self-confidence, improved chewing or biting, improved speech, and improved ability to comfortably close lips.

Indications:

The Invisalign System is indicated for the orthodontic treatment of malocclusion. The optional mandibular advancement feature(s) positions the patient's jaw to address skeletal malocclusion.

After initial distribution in the dental professional's office, the Invisalign System is intended for home use.

Contraindications:

The Invisalign system is contraindicated for use in patients with active periodontal disease.

Note: this does not preclude, at the discretion of the doctor, treatment of patients with treated, controlled, or arrested periodontal disease and free of active disease.

Warnings

1. In rare instances, some patients may be allergic to the plastic aligner material (e.g., plastic, coating material) including aligners

with occlusal blocks material.

2. In rare instances, patients with hereditary angioedema (HAE), a genetic disorder, may experience rapid local swelling of subcutaneous tissues including the larynx. HAE may be triggered by mild stimuli including dental procedures.
3. Orthodontic appliances, or parts thereof, may be accidentally swallowed or aspirated and may be harmful. **Note:** If any of the above should occur, patient should discontinue use, seek medical attention immediately as necessary, and notify the Invisalign treating doctor. The doctor is required to notify Align Technology and the patient may also notify Align Technology, if needed.
4. Orthodontic treatment, including clear aligner treatments, is not effective in the movement of dental implants and should not be used to move implants. Treatment plans that involve implant movement may result in an implant failure.
5. Invisalign treatment, including aligners with mandibular advancement features, is not tested nor intended for treatment of respiratory disorders, sleep apnea, bruxism, and/or temporo-mandibular disorders (TMD).
6. Erupting, exfoliating and/or missing teeth under the mandibular advancement features may affect the integrity of the feature. In addition, any erupting teeth which are located under the features may not complete eruption during treatment with the features.

Precautions:

Treatment considerations

1. Aligner treatment alone may not be adequate to achieve the desired result in patients with certain orthodontic issues. Your doctor will determine whether supplemental treatment is necessary and will discuss any associated risks and healing requirements prior to starting aligner treatment.
2. In patients with small, unusually - shaped or missing teeth, aligner retention and treatment duration may be impacted.
3. Dental restorations (which replace or restore part or all of a patient's tooth – for example, crowns or bridges) may become dislodged and require re-cementation or, in some instances, replacement.
4. A tooth that has been injured (for example, from a fall or blow) may be aggravated during orthodontic treatment. In rare instances, the life of the tooth may be reduced, the tooth may require additional dental treatment such as endodontic (root canal) and/or additional restorative work, and/or the tooth may be lost.
5. The root lengths of the teeth may be shortened (this is known as 'root resorption') during orthodontic treatment, causing a threat to longevity of the teeth.
6. Do not alter your aligners. The aligners should cover all your teeth; otherwise, 'supra-eruption' (in which an uncovered tooth can move above the 'occlusal' or chewing surface of the other teeth) may occur.
7. In rare instances, problems in the temporo-mandibular or jaw joint (temporo-mandibular disorder or dysfunction (TMD)) may result in joint pain, headaches, or ear problems. For patients in treatment with mandibular advancement features (which are used to shift the lower jaw forward), problems in the jaw joint may be exacerbated.
Note: In patients with limited vertical opening or active TMD, doctor may prefer to prescribe precision wings instead of occlusal blocks.
8. Certain medical conditions and use of certain medications may affect orthodontic tooth movement and treatment outcomes.
9. Do not cut, remove, or try to shave off the blocks from the aligners. This will impact the product integrity and prevent the aligners with mandibular advancement features such as occlusal blocks work as intended.

Attachments and interproximal reduction (IPR)

1. Attachments (bumps which are attached to your teeth during aligner treatment to help your aligners move your teeth) may be temporarily bonded to one or more teeth during the course of treatment to help hold your aligners in place and to aid tooth movement.
2. Attachments may fall off and require replacement. Failure to have attachments replaced may impact a patient's treatment duration or treatment outcome.
3. All attachments should be removed when your aligner treatment is complete. Failure to remove attachments at the end of treatment could lead to tooth decay (cavities) and soft tissue irritation.
4. Interproximal reduction (IPR) (filing in between your teeth) may be prescribed to create space as part of treatment. You may experience temporary sensitivity and your doctor will discuss any associated risks and remedies.

Treatment experience

1. Aligners with cracks the size of the height of the tooth should not be used.
2. Dental tenderness or sensitivity may be experienced following initial aligner placement and after switching to each new aligner in the series.
3. Patients may experience a temporary loosening of their teeth during the treatment.
4. The product may temporarily affect speech and may result in a lisp, although any speech impediment associated with clear

aligner treatment usually disappears within one or two weeks.

5. A temporary increase in salivation or dryness of mouth may occur.
6. Gums, cheeks, or lips may be scratched or irritated by the product and its associated features.
7. Tooth decay (cavities), periodontal disease (gum disease), and permanent markings from stains and decalcification (white spots) may occur if patients do not brush and floss their teeth properly during treatment or if they consume foods or beverages containing sugar or acids while wearing aligners.
8. The bite may change during treatment; this may result in temporary discomfort.
9. At the end of treatment, the bite may require adjustment by the doctor.
10. Teeth may shift position (relapse) after treatment. Consistent wearing of retainers at the end of treatment should reduce this tendency.
11. For aligners with mandibular advancement feature, biting directly or grinding on the features or positioning them incorrectly may increase the chance of aligner deformation or breakage, which may lengthen the treatment.
12. Prolonged biting or grinding on and/or improper positioning of occlusal blocks may also lead to soreness in the muscles or jaw joint.

Aligner care and compliance

1. Aligners are intended to be removed only for eating, drinking, and oral care. Unnecessary and excessive removal and insertion of aligners may prematurely damage and/or break the aligners.
2. Failure to properly store, insert or remove aligners can lead to damage. Patients should read all instructions carefully and should use proper techniques for aligner insertion and removal. Aligners should always be stored in a cool, dry place and in their applicable cases. Attempting to store aligners with mandibular advancement features in a standard Invisalign Aligner Case may damage the aligners.
3. Keep aligners away from hot water and harsh chemicals and carefully follow the patient instructions below. Do NOT use denture cleaners to clean aligners with mandibular advancement features or soak them in mouthwash. Do NOT use Invisalign Aligner Cleaning Foam, Invisalign Aligner Cleaning Spray, or Invisalign UV Sanitizing Case to clean aligners with mandibular advancement features. These products may damage the surface of the aligner, causing it to become dull, more visible, and may affect bonding of blocks to aligners.
4. Keep aligners out of reach of young children and pets. For children in treatment, a supervising adult should monitor aligner use, storage and care.
5. Tobacco use, smoking, and/or vaping while wearing aligners may lead to discoloration or damage. In general, tobacco use has been shown to increase gum disease and to delay tooth movement during orthodontic treatment.
6. Orthodontic treatment (including aligner treatment) may impair the health of the bone and gums which support the teeth and may aggravate the gums.
7. The length and success of treatment depends on the patient's compliance, cooperation in keeping appointments, maintaining good oral hygiene, avoiding loose or broken appliances and following the doctor's instructions carefully.
8. Failure to wear the appliance for the prescribed number of hours per day and/or not using the product as directed by the doctor can lengthen the treatment time and affect the ability to achieve the desired results. The Compliance Indicator feature can be order to approximate gauge of patient wear time by the doctor.
9. Regardless of the patient's compliance, due to the variation in the size and shape of the teeth and other factors, achievement of an ideal result may not always be possible. Restorative dental treatment may be indicated.

Wearing and using your Invisalign® aligners

Here are some tips to help ensure proper use and avoid damaging your aligners.

Always remember to

1. Wear your aligners per your doctor's instructions, usually 20-22 hours per day.
2. Wash your hands thoroughly with soap and water before handling your aligners.
3. Handle only ONE aligner at a time.
4. Rinse your aligners when removing them from the packaging.

Use proper technique described below when inserting and removing aligners to avoid damaging aligners.

Aligner insertion

1. Prior to each insertion and before bed, inspect aligners for cracks and deformation. Severe cracks or deformities should be reported to the doctor immediately. Caution on aligners with cracks the size of the height of the tooth, or occlusal blocks that

are deformed or that appear to be debonding from the aligner (figure 1) should not be used.

Note: Parent or guardian should inspect aligners of pediatric patients.

2. Make sure you have the proper aligner—the upper for your top teeth and the lower for your bottom teeth.
3. To help avoid confusion, each aligner is engraved with your unique case number, a “U” for upper and an “L” for lower, followed by the stage number (figure 2).
4. You may insert either the upper or lower aligner first. When inserting each aligner, gently push the aligners over your front teeth. Then apply equal pressure, using your fingertips, to the tops of your left and right molars (back teeth) until the aligner snaps into place (figure 3).
5. Do NOT bite your aligners into position. This may damage them. For aligners with precision wings (figure 4 and 5) do NOT bite down or grind on precision wings or occlusal blocks (figure 6 and 7). The lower precision wings should be positioned in front of the upper precision wings when the aligners are in the mouth; with occlusal blocks, the lower blocks should be positioned in front of the upper blocks.

Aligner removal

1. On one side of your mouth, use your fingertip on the inside of your back molar to slowly pull the aligner from your molars (figure 8).
2. Repeat this process on the other side of your mouth before trying to completely remove the aligner. (figure 9)
3. Once aligner is disengaged from the molars on both sides of your mouth, you should be able to slowly work your way forward gently prying the aligner away from your teeth with your fingertips. (figure 10)

Note:

- Immediately rinse aligner with water, shake off excess water, and store your aligners in the protective case provided with your starter kit.
- Aligners are intended to be removed only for eating, drinking, smoking, and oral care. Unnecessary and excessive removal/insertion of aligners with mandibular advancement features may prematurely damage and/or break the aligners.
- Take care in removing your aligners, especially if multiple attachments are being used.
- Do not use excessive force to bend or twist an aligner to get it off.
- DO NOT use any sharp object to remove your aligners. Consider an aligner removal tool if needed.
- Consult with your doctor if your aligners are extremely difficult to remove.

Daily care and maintenance of your invisalign aligners

1. Clean your aligners prior to each insertion. Use a soft bristle toothbrush with water and a small amount of toothpaste. You may find it easiest to clean the outside of your aligners by brushing them while they are still on your teeth (figure 11), then remove your aligners to clean the inside surfaces (figure 12).

Note: Be sure to rinse each aligner thoroughly with water after each cleaning. We also recommend that you use the Invisalign Cleaning Crystals once a week, or as needed.

2. Prior to each insertion and before bed, inspect aligners for cracks and deformation. Severe cracks, deformities or block debonding should be reported to the doctor immediately. Aligners with cracks the size of the height of the tooth or occlusal blocks that are deformed, or that appear to be debonding from the aligner should not be used.

Note: Parent or guardian should inspect aligners of pediatric patients.

3. Do NOT use denture cleaners to clean aligners including aligners with mandibular advancement features or soak them in mouthwash. Do NOT use Invisalign Aligner Cleaning Foam, Invisalign Aligner Cleaning Spray, or Invisalign UV Sanitizing Case to clean aligners with mandibular advancement features. These products may damage the surface of the aligner, causing it to become dull and more visible and may affect bonding of blocks to aligners.
4. Invisalign Aligner Cases may be hand washed with soap and water. Do not use dish washing machines to clean the cases.

Proper oral hygiene

1. Remove your aligners for eating and drinking. (You do not need to remove your aligners to drink cool water.)
2. Brush and floss your teeth after each meal or snack prior to re-inserting your aligners. If you don't have access to your cleaning system or a toothbrush, you can simply rinse your mouth, and then clean your aligners by holding them under warm running water. It's not the best way to clean, but it works in a pinch. Be sure to thoroughly clean your aligners at your earliest convenience.
3. If you have any questions regarding hygiene techniques, please consult your doctor.
4. Regular dental checkups and cleaning are recommended for the continued health of your teeth and gums.

Storing your Invisalign aligners

1. Invisalign Aligner Case(s) are included in your starter kit. We recommend that you store your Invisalign aligners in this case when they are not in your mouth (figures 13 and 14). This will protect them from loss and damage. Keep aligner cases clean by hand washing them with soap and water. Do not wash them in dish washing machines.
2. The case represented on Figure 14 is designed for aligners that have "wings" on the sides or "blocks" on the chewing surface of aligners. Trying to place these aligners into the standard Invisalign Aligner Case shown in Figure 13 may damage the aligners, which could impact treatment efficacy

Always keep the most recently used aligners (or more stages as directed by your doctor). If your current aligner is lost or broken, your doctor may recommend that you temporarily go back one stage while a replacement is being made.

DEVICE DISPOSAL

After the treatment is complete as deemed by the doctor, Invisalign aligners should be disposed of in compliance with national and local laws.

Note: Keep all of your older aligners in a clean plastic bag, or as instructed by your doctor. Your doctor may recommend that you temporarily go back one stage and use one of the older aligners while a replacement is being made. Keep them out of reach of small children or pets.

Frequently asked questions and answers

For patients undergoing mandibular advancement, why do some of my aligners have 'wings' or 'blocks' on them?

The 'wings' on the sides or the 'blocks' on the chewing surface of aligners are called precision wings or occlusal blocks, and are especially designed to position the lower jaw forward to help correct the bite.

What should I do if my new Invisalign aligner doesn't snap onto my teeth?

Minor discrepancies between the new aligner and the current tooth position are normal, since the teeth need time to conform to the new aligner position. In the event of significant problems with aligner fit, inform your doctor.

Is it OK to drink hot or cold beverages while wearing Invisalign aligners?

Except for cool water, we recommend that you do not drink while wearing the aligners. This is to avoid formation of cavities and stains, or warping of the aligners with hot drinks and hot water.

What if I lose or break an attachment?

In the event that an attachment is lost or broken, you should immediately contact your doctor.

Note: If you have additional questions after reading this pamphlet, please consult with your doctor.

Customer Information/User Information

In case of a serious incident that directly or indirectly led, might have led or might lead to any of the following: a) the death of a patient, user or other person, b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health, c) a serious public health threat please immediately without any delay inform Align Technology and your local competent authority.

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Consult instructions for use or consult electronic instructions for use

[illegible]

Use-by date

• تاريخ انتهاء الصلاحية • Срок на годность • Datum isteka roka valjanosti • Data de caducitat • Datum spotřebny • Sidste anvendelsesdato • Houdbaarheidsdatum • Kõlblikkuskuupäev • Viimeinen käyttöpäivä • Date limite d'utilisation • Verfallsdatum • Ημερομηνία λήξης • Felhasználhatósági időpont • Data di scadenza • Deriguma terménsh • Tinkamumo naudoti terminas • Utlopsdatum • Termin ważności • Utilitar de preferență ate • Data limită de utilizare • Годен до • Rok upotrebe • Spotřebuje do • Rok uporabnosti • Fecha de caducidad • Förbrukningsdatum • Son kullanna tarihi • Термін придатності

Unique device identifier

مفرد الجوار الفريد – Unikálnen identifikátor na izdelatelo / Jedinstvena identifikacija proizvođa / Identifikator único del
 الجوار – Jedinečný identifikátor prostředku / Unik ehensidifikator / Unieke apparaatidentificatie / Seadme kordumatu
 identifikatsioon / Yksilöllinen laitteistunus / Identifiant unique du dispositif / Eindeutige Gerätekennung / Монодік
 αναγνωριστικό / μονοσήμαντικό προϊόντος / מזהה יחיד / Egyedi eszközazonosító / Identificativo unico del dispositivo
 / Unikals ierices identifikators / Unikals prietaiso identifikatorius / Unik ehensidifikator / Niepowtarzalny kod identyfikacyjny
 wyrobu (ang. unique device identifier, UDI) / Identificador único do dispositivo / Identifikator único do dispozitívu / Унікальний
 ідентифікатор виробу / Jedinstveni identifikator sredstva / Jedinstven identifikator / Identifikator único do produto /
 pripomočka / Identificador único del producto / Unik ehensidifikatore / Benzersiz cihaz tanımlayıcı / Унікальний ідентифікатор
 пристрою

Country of manufacture combined with the date of manufacture

3. в. الجمع مع ذكر تاريخ التصنيع : Комбинация от др-жава на производство и дата на производство • Zemlja proizvodnje i datum proizvodnje • Pays de fabrication combinat avec la date de fabrication • Země a datum výroby • Produktionslandet kombinert med produktionsdatoen • Land van vervaardiging gecombineerd met de productiedatum • Tootariik koos tootmise kuupäevaga • Valmistusmaa ja valmistuspäivä • Pays de fabrication combiné à la date de fabrication • Herstellingsland zusammen mit dem Herstellungsdatum • Γαράτση κατασκευής σε συνδυασμό με την ημερομηνία κατασκευής • Gyártó ország a gyártási dátummal kombinálva • Paese di fabbricazione e data di fabbricazione • Ražotājvalsts un ražošanas datums • Pagamino šalis kartu su pagamino data • Produktionsland sammen med produktionsdatoen • Kraj produkcji wraz z datą produkcji • Pays de fabrico combinado com a data de fabrico • Tara de fabricație combinată cu data de fabricație • Страна, в которой произведен продукт, и дата изготовления • Zemlja proizvodnje u kombinaciji sa datumom proizvodnje • Країна výroby в комбінації з датумом výroby • Država izdavele in datum proizvodnje • El país de fabricación de los productos y la fecha de fabricación • Tillverkningsland kombinert med tillverkningsdatum • Üretim tarihi ile birlikte üretim ülkesi • Країна виробництва в поєднанні з датою виробництва

Patient number

1. **Номер на пациента** • Broj pacijenta • Número de pacient • Číslo pacienta • Patientnummer • Páciéntsám • Patsiendi number • Potilasnumero • Numéro du patient • Patientennummer • Αριθμός ασθενούς • מספר המטופל • Betegszám • Numero del paciente • Pacienta numurs • Paciento numeris • Patientnummer • Numer pacijenta • Número do paciente • Numărul pacientului • Номер пациента • Broj pacijenta • Číslo pacienta • Stevilka bolnika • Número de paciente • Patientnummer • Hasta numarasi • Номер пациента

	Single patient multiple use استخدام متعدد لمريض واحد • Един пациент, многократна употреба • Jedan pacijent, višekratna upotreba • Ūs mŭltiple en un sol pacient • Jeden pacient, viacnásobné pouŭitie • Enkelt patient, flere anvendelser • Enkele patiŕnt, meervoudig gebruik • Ūks patiŕnt, mitu kasutuskorda • Yksi potilas, useita kŕyttŕkertoja • Patient unique, utilisation multiple • Mehrfachanwendung fŭr einen einzelnen Patienten • Ένας ασθενής, πολλαπλή χρŕση • חֶדֶד לִשְׁמוֹשׁ בְּפַעַם בְּסוּפֵלוֹ • Egyetlen beteg, tŕbbszŕri hasznŕlat • Paziente singolo, piŭ utilizzi • Viens patients, vairŕkŕi lŕetajumi • Vienas pacientas, daugkartinis naudojimas • Enkeltpatient, flerbruk • Do wielokrotnego uŭcinka u jednego pacjenta • Um ŭnico paciente, vŕrias utilizaçŕes • Patient unic, utilizare multiplŕ • Многократное использование для одного пациента • Za višekratnu upotrebu od strane jednog pacjenta • Jeden pacient, viacnŕsobnŕ pouŭitie • Za enega bolnika in za veċkratno uporabo • Un paciente, varios usos • En patient, flera anvŕndningar • Tek hasta, birden fazla kullanim • Один пациент, багатаразове використання
EC REP	Authorized representative in the European Community/European Union الممثل المفوض في الجماعة الأوروبية/الاتحاد الأوروبي • Упълномощен представител в Европейската Общност/Европейския съюз • Ovlaŕteni zastupnik za Europsku zajednicu / Europsku uniju • Representant autorizat a la Comunitat Europea/Uniŕo Europea • Zplnomocnenŕy zŕstupce v Evropskŕm spoŕlenstvŭ / Evropskej unii • Autoriseret reprŕsentant i Det Europŕiske Fŕllesskab/Den Europŕiske Union • Gemachtigd vertegenwoordiger in de Europese Gemeenschap/Europese Unie • Volitatus esindaja Euroopa Ūhenduses / Euroopa Liidus • Valtuutettu edustaja Euroopan yhteisŕssŕ / Euroopan unionissa • Reprŕsentant autorisŕ pour la Communautŕ europŕenne/l'Union europŕenne • Autorisierte Vertretung in der Europŕischen Gemeinschaft/Europŕischen Union • Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα/Ευρωπαϊκή Ένωση • נציג מורשה בקהילה האירופית/באיחוד האירופי • Meghatalmazott képviselŕ az Eurŕpai Kŕzŕssŕgben/Eurŕpai Uniŕban • Rappresentante autorizzato nella Comunitŕ europea/Unione europea • Pilnvarotais pŕrstŕvis Eiropas Kopienŕ/Eiropas Savienibŕ • Įgaliojatis atstovas Europos Bendrijŕje / Europos Sąjungaŕje • Autoriseret representant i Det Europŕiske fŕllesskab / EU • Upowaŭniony przedstawiciel we Wspŕlnocie Europejskiej / Unii Europejskiej • Representante autorizado na Comunidade Europeia/Uniŕo Europeia • Reprezentant autorizat in Comunitatea Europeanŕ/Uniunea Europeanŕ • Уполномоченный представитель в Европейском сообществе / Европейском союзе • Ovlaŕŕeni predstavnik u Evropskoj zajednici/Evropskoj uniji • Zplnomocnenŕy zŕstupca v Eurŕpskom spoŕlenstve/Eurŕpskej ŭnii • Pooblaŕŕŕeni zastopnik v Evropski skupnosti/Evropski uniji • Representante autorizado en la Comunidad Europea/Uniŕn Europea • Auktoriserad representant i Europeiska gemenskapen/Europeiska unionen • Avrupa Topluŭuĝu'nda/Avrupa Birliĝi'nde yetkili temsilci • Уповноважений представник у Европейському Спŕтовариствŭ / Европейському Союзŭ
CH REP	Swiss authorized representative الممثل المسموحي السويسري • Упълномощен представител за Швейцария • Ovlaŕteni zastupnik za Ŗvicarsku • Representant autorizat suis • Zplnomocnenŕy zŕstupce ve Ŗvicarsku • Autoriseret reprŕsentant i Schweiz • Gemachtigde vertegenwoordiger van Zwitserland • Sveitsi volitatus esindaja • Valtuutettu edustaja Sveitsissŕ • Reprŕsentant autorisŕ pour la Suisse • Autoriserter Vertreter in der Schweiz • Εξουσιοδοτημένος αντιπρόσωπος στην Ελβετία • נציג מורשה בשווייץ • Svŕjci meghatalmazott képviselŕ • Rappresentante autorizzato in Svizzera • Pilnvarotais pŕrstŕvis Svŕicŕ • Įgaliojatis atstovas Ŗveicŕje • Įgaliojatis atstovas Ŗveicŕje • Sveitsi auktoriseret representant • Upowaŭniony przedstawiciel w Ŗwajcarii • Representante autorizado na Suia • Representant autorizat in Elvetia • Уполномоченный представитель в Швейцарии • Ovlaŕŕeni predstavnik u Ŗvajcarskoj • Pooblaŕŕŕeni zastopnik v Ŗvajcarski skupnosti • Pooblaŕŕŕeni zastopnik v Ŗvici • Representante autorizado suizo • Auktoriserad representant i Schweiz • Isvŕŕŕe yetkili temsilcisi • Уповноважений представник у Швейцарii
	Importer المسموحي • Вноситель • Uvoznik • Importador • Dovoze • Importŕ • Importeur • Importtija • Maahantuoja • Importateur • Importeur • Εισαγωγέας • יבואן • Importŕ • Importatore • Importŕtŕjs • Importuotojas • Importŕ • Importer • Importador • Importator • Импортѕр • Uvoznik • Dovoza • Uvoznik • Importador • Importŕ • lthalatŕ • Импортѕр
	Знак вŕдповŕдностŭ Технŕчному регламенту

Medical devices-Symbols to be used with medical device labels, labelling and information to be supplied-Part 1: General requirements.

For professional dental use – **Rx only**

align™

Align Technology, Inc.
2820 Orchard Parkway
San Jose, CA 95134, USA

Align Technology Switzerland GmbH,
Suurstoffi 22
6343 Risch-Rotkreuz, Switzerland

Align Technology B.V.
Herikerbergweg 312
1101 CT Amsterdam, The Netherlands

EU Importer

Align Technology (Poland) Sp. z o.o.
Bierutowska 57-59
51-317 Wrocław, Poland

UK responsible person:
Align Technology UK Ltd.
2800 The Crescent, Solihull Parkway
Birmingham Business Park, Solihull
B37 7YL, England, United Kingdom

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Patent www.aligntech.com/patents

© 2024 Align Technology, Inc. Align, Invisalign, the starflower mark, and invis are trademarks of Align Technology, Inc.
224244 Rev B (2024-05)

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