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

# Firmware update of the Axon-Bus Prosthetic System with AxonSoft


Ensure that the following conditions are fulfilled before the firmware update is carried out:

- Completely assembled Axon-Bus Prosthetic System with Axon-Bus gripping unit
- Fully charged battery of the Axon-Bus Prosthetic System
- Installed AxonSoft software on the computer to be used
- BionicLink PC 60X5 connected on the computer to be used and all corresponding drivers are installed
- Axon-Bus Prosthetic System and BionicLink PC 60X5 are within close proximity to one another (within the Bluetooth range)
- Downloaded Firmware update

1. Install the downloaded firmware update and follow the installation wizard instructions.
2. Start AxonSoft in the Ottobock Data Station.
3. Switch off the Axon-Bus Prosthetic System.
4. Press the button in the charging receptacle of the Axon-Bus Prosthetic System for approximately 4 seconds to switch into the Bluetooth connection mode.
  - The Light Emitting Diode (LED) in the charging receptacle lights up blue.
5. In the **Connection** tab, click the **Connect** button to establish a connection to the Axon-Bus Prosthetic System.
6. Select the desired Axon-Bus Prosthetic System.
7. Click the **OK** button for confirmation.
  - The LED in the charging receptacle flashes blue.
  - The following window appears before the update starts:

8. Click the **Yes** button to install the update.
  - The Update process starts.
  - **Note: Do not touch or disturb the Axon-Bus Prosthetic System and/or the BionicLink PC 60X5 while updating, as this could terminate the connection.**
9. The following information will be displayed after the update process is successfully completed:

- **Note: If an error message appears, recheck the aforementioned firmware update conditions (see begin of this chapter) and repeat the updating process.**
10. If the error message still persists, please contact Reinhard Walter at: [reinhard.walter@ottobock.com](mailto:reinhard.walter@ottobock.com).