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

After an hour of stimulation, patients were already able to move more freely and safely. Four weeks of regular use helped them maintain this rapid improvement.

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Patients with pain reported that their symptoms were rapidly reduced by an hour in the Suit. After four weeks of regular, near full-body stimulation, they continued to feel significantly less pain.

Exopulse Registry

Open-label study of group-level response in 38 patients with CP (adult + pediatric), MS, or stroke.

- All patients: GMFCS 1-3, BBS score < 45
- CP: n=15 (9 adult, 6 pediatric), mean age 20.9
- MS: n=15, mean age 50.3
- Stroke: n=8, mean age 43.6

Key features.

An innovative design tailored to your patient.

Technology

- **Near full-body stimulation:** The unique wearable configuration is designed to manage spasticity and related pain in both the upper and lower limbs.
- **58 embedded electrodes:** Stimulation through each channel can be customized for the patient, enabling muscle-by-muscle targeting of up to 40 muscle groups.
- **Convenient control unit:** The Suit's activation and programming hub is easily accessible for both clinicians and patients.

Garments

- **Jacket and pants:** Suit is made up of two wearable pieces that are easy for patients to put on at home.
- **Construction:** Both garments combine comfortable tailoring with durable, reinforced design, optimized for regular use.
- **Material:** The jacket and pants are made of a comfortable, breathable synthetic blend.
- **Machine washable:** Both garments (without the control unit) can be cleaned up to 25 times in a home washer, or regularly in a steam cabinet.

Using the Suit.

Convenient, at-home spasticity management.

Here's a step-by-step summary of how clinics like yours can make the Suit part of your patients' care.

- 1 Get training for your team.** Before you start fitting patients, Ottobock experts will guide you through setting up the components, programming the Suit, and coaching your patients on how to maintain it.
- 2 Get your patient started.** First, your patient will visit your clinic so you can show them how to wear and take care of their Suit. During that visit, you'll also customize the stimulation pattern for the patient's unique needs.
- 3 Send them home with their Suit.** Once your patient and/or their caregivers have been trained on the Suit, they can take it home to use as directed. Most spasticity patients should wear it every other day.
- 4 Monitor their progress.** Check in regularly to ensure patients and caregivers are using their Suit properly and getting the results they want to see. If needed, you can make more adjustments to their stimulation.
- 5 Take their standard treatments further.** Regular use of the Suit may help your patient get more out of their ongoing physical therapy or exercise regimen. Unless contraindicated, patients can also continue using prescribed treatments like oral medications or BoTN injections (as directed by their physician).

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*Based on change in Berg Balance Scale (BBS) score in 38 patients with impaired balance and an increased fall risk (n=15/15/8 CP/MS/stroke, baseline BBS score < 45). Patients in all cohorts reported significant improvement in BBS score after sixty minutes of stimulation (T1) and after 4 weeks of stimulation every other day (T2).

**Based on change in EQ5D pain levels in a subset of 32 patients who reported spasticity-related pain at baseline (n=10/15/7 CP/MS/stroke, baseline EQ5D pain > 1). Patients in all cohorts reported significant improvement in EQ5D pain score after sixty minutes of stimulation (T1) and after 4 weeks of stimulation every other day (T2).

References: 1. Elbasiouny SM, Moroz D, Bakr MM, Mushahwar VK. Management of Spasticity After Spinal Cord Injury: Current Techniques and Future Directions. *Neurorehabil Neural Repair*. 2010; 24(1):23-33. DOI: 10.1177/1545968309343213. 2. Rabchevsky AG, Kitzman PH. Latest approaches for the treatment of spasticity and autonomic dysreflexia in chronic spinal cord injury. *Neurotherapeutics*. 2011; 8(2):274-282. DOI: 10.1007/s13311-011-0025-5. 3. Stein C, Fritsch CG, et al. Effects of Electrical Stimulation in Spastic Muscles After Stroke: Systematic Review and Meta-Analysis of Randomized Controlled Trials. *Stroke*. 2015; 46(8):2197-2205. DOI: <https://doi.org/10.1161/STROKEAHA.115.009633>. 4. Bosques G, Martin R, et al. Does therapeutic electrical stimulation improve function in children with disabilities? A comprehensive literature review. *J Ped Rehab Med*. 2016; 9(2):83-99. DOI: 10.3233/PRM-160375. 5. Etoom M, Khraiweh Y, et al. Effectiveness of Physiotherapy Interventions on Spasticity in People With Multiple Sclerosis: A Systematic Review and Meta-Analysis. *Am J Phys Med Rehab*. 2018; 97(11):793-807. DOI: 10.1097/PHM.0000000000000970. 6. Exopulse Registry Clinical Research Report; Data on File.

© 2022 Ottobock. All rights reserved.

Notification: The information contained in this document is only intended for an audience outside the United States. Please note that the products described herein are not FDA approved.

Not all products and services are registered or available for sale in all countries.

Ottobock SE & Co. KGaA
Max-Näder-Straße 15
37115 Duderstadt
Germany
www.ottobock.com