ottobock.

Otkrij svoje supermoći

Oslobodite se spastičnosti i
bola vezanog za spastičnost
pomoću *Exopulse Mollii Suit*.

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Skenirajte da biste pogledali Louisinu priču.

”

The Exopulse Mollii Suit empowers me to relax my muscles and get lost in the music.

Louisa,
Multiple sclerosis

More information about the *Exopulse Mollii Suit*.

Have more questions or want to talk to an expert?
Find an Ottobock clinical partner near you at ottobock.com.

Where can I or my loved one try out the Suit?

The Suit is currently available at select clinics and rehabilitation facilities. The certified clinicians at these locations are trained on using the Suit to manage spasticity in patients with conditions like CP, MS, and stroke. To find an Ottobock clinical partner near you, visit ottobock.com.

Do I need a prescription to use the Suit?

Depending on where you're located, or which clinic provides your Suit, you may or may not need a prescription to use it. If you're interested in trying out the *Exopulse Mollii Suit*, ask your doctor if you need a prescription to do so (or to seek reimbursement for using the Suit).

Does the Suit have any side effects?

Patients wearing the Suit have rarely experienced some mild, temporary side effects of the stimulation they receive. These side effects typical of electrical stimulation include minor skin irritation, surface-level skin redness, and mild burns.

Can I keep using other treatments?

In most cases, yes, but always consult your doctor to ensure that this is appropriate for you or your family member. With your doctor's approval,

the *Exopulse Mollii Suit* should ideally be used in combination with physiotherapy, training, or physical activity. Discuss the Suit with your doctor if you or your loved one are currently receiving botulinum toxin injections or have been prescribed oral medication for spasticity. Do not use the Suit if you or your family member has an implanted medical device.

How much does the Suit cost? Is it covered by insurance or my national health system?

The cost of the Suit may vary depending on your age, your diagnosis, and other factors. Financial coverage or reimbursement for using the Suit will differ depending on where you're located. For more information about the potential cost of the Suit, contact an Ottobock clinical partner near you.

How do I take care of the Suit? Is it machine-washable?

The jacket and pants (not the control unit) are machine-washable up to 25 times at 40°C/104°F. It may also be cleaned in a steam cabinet as often as necessary. The control unit may NOT be placed in a washing machine or steam cabinet. Doing so will void the two-year warranty.

Has the Suit been evaluated in any clinical studies?

The Suit's effectiveness has been evaluated in several different clinical research programs, including both independent and Ottobock-sponsored studies. If you'd like to learn more about the results of these studies, contact an Ottobock clinical partner near you.

Caution

Do not use the *Exopulse Mollii Suit*:

- If you have implanted electronic medical devices or equipment that can be disrupted by magnets (e.g., shunts)
- Together with electronic life-support equipment or high-frequency operation equipment
- Together with ECG-equipment

The *Exopulse Mollii Suit* may potentially disrupt the function of the above-mentioned types of equipment. All incorrect usage is at the user's own risk.

Stimulation should not be applied:

- Over swollen, infected, or inflamed areas or skin eruptions (e.g., phlebitis, thrombophlebitis, varicose veins, etc.)
- Over the neck or mouth. Severe spasms of the laryngeal and pharyngeal muscles may occur, and the contractions may be strong enough to close the airway or cause difficulty in breathing
- Transthoracically, in that the introduction of electrical current into the heart may cause cardiac arrhythmias
- Transcerebrally

Do not use the *Exopulse Mollii Suit* without consulting a doctor in connection with:

- Epilepsy
- Cardiovascular diseases
- Malignancy (cancer)
- Infectious diseases
- Fever
- Pregnancy
- Skin disease, rashes, or other skin problems
- Another medical device or medical treatment

Using the *Exopulse Mollii Suit* with the above can expose the user to unnecessary risk and is therefore at the user's own risk. Always consult a certified *Exopulse Mollii Suit* clinician or your doctor.

Care and maintenance

Machine-washable

The jacket and pants may be cleaned up to 25 times in a home washer, or regularly in a steam cabinet.

Warm water only

Do not clean the jacket and pants in cold or hot water.

Do not wash the control unit

This may void the Suit's two-year warranty.

Javite se da isprobate.

Spremni za probu? Proba odela je trenutno dostupna u prostorijama Otto Bock Sava d.o.o. u Kragujevcu, u Industrijskoj bb.

Sertifikovani Mollii eksperti u Ottobock Srbija su posebno obučeni za korišćenje odela za upravljanje spastičnosti kod pacijenata sa stanjima kao što su CP, MS i moždani udar. Da biste videli da li je odelo pravo za vas, kontaktirajte Ottobock u vašoj oblasti i zatražite probu odela.

Obaveštenje: Informacije sadržane u ovom dokumentu namenjene su samo publici van Sjedinjenih Država. Ovaj proizvod nije naveden ili odobren od strane FDA. Nisu svi proizvodi i usluge registrovani ili dostupni za prodaju u svim zemljama.

Ottobock SE & Co. KGaA
Max-Näder-Straße 15
37115 Duderstadt
Germany
www.ottobock.com