RECOMMEND PRODUCTS FOR YOUR PATIENTS, FAMILY AND FRIENDS THROUGH THE



Electric Rechargeable Systems



TRANSFORMATIONAL GUM HEALTH SYSTEM

- iO[™] Electric Toothbrush
- PRO-HEALTH™ Gum Detoxify™ Toothpaste
- PRO-HEALTH™ Multi-Protection Clean Mint Mouthwash
- PRO-HEALTH $^{\text{\tiny{M}}}$ Deep CleanFloss

Brush head refills include:





Gentle Care



GINGIVITIS SYSTEM

- GENIUS™ X Professional Exclusive Electric Toothbrush
- PRO-HEALTH[™] Gum Detoxify[™] Toothpaste
- PRO-HEALTH™ Multi-Protection Clean Mint Mouthwash
- PRO-HEALTH[™] Deep Clean Floss

Brush head refills include:



 $\mathsf{CrossAction}^{\scriptscriptstyle\mathsf{TM}}$



FlossAction™



Sensitive Gum Care



Interproximal Clean



HEALTHY WHITE SYSTEM

- PRO 1500 Electric Toothbrush
- 3D WHITE™ Brilliance Toothpaste
- 3D WHITE™ Glamorous White Mouthwash
- PRO-HEALTH™ Advanced Floss

Brush head refills include:



 $\mathsf{CrossAction}^{\scriptscriptstyle\mathsf{TM}}$



RECOMMEND PRODUCTS FOR YOUR PATIENTS, FAMILY AND FRIENDS THROUGH THE



Crest® Whitening

\$179.97

per case

of 3



CREST 3D WHITESTRIPS + LED LIGHT

1 hour for visible results.

Professional-level teeth

whitening in just 7 treatments.

Dental-Professional Exclusive

WHITENS TEETH 1.5% BETTER

LIST BETTER

ENING ACCELERATOR

CREST WHITENING EMULSIONS + LED LIGHT

Noticeably whiter teeth after 14 treatments.

Highly active micro-droplets with over 20% hydrogen peroxide suspended in a hydrating base.

3% peroxide overall.

30 minutes + LED light time. Apply up to 4 times a day. Do not rinse or brush for 30 minutes after use. \$180.00 PETERNING MULSIONS POR CASE OF 8

CREST WHITENING EMULSIONS ON-THE-GO

For touch-ups to maintain a white smile with the convenience of on-the-go use.

Highly active micro-droplets with over 20% hydrogen peroxide suspended in a hydrating base.

3% peroxide overall.

I minute (Apply & Go). Apply up to 4 times a day. Do not rinse or brush for 30 minutes after use.

MODE OF ACTION & DELIVERY SYSTEM

EXPECTED RESULTS

11% hydrogen peroxide.

APPLICATION TIME

1 hour + LED light time. Maximum once a day.

- All Crest at-home whitening products are enamel safe.
- Results depend on different factors, including a patient's age, how yellow their teeth are, and how often they apply the product.

 Most patients see teeth whitening results after just a few uses.
- Duration of results varies, depending on the patient's use of teeth-discoloring beverages like wine, coffee, tobacco, soda, etc.
- Remember, at-home whitening products won't work on caps, crowns, veneers, partials, fillings, laminates, or other similar dental work.

Crest + Oral: