

RECOMMEND PRODUCTS FOR YOUR PATIENTS, FAMILY AND FRIENDS THROUGH THE

Crest + Oral-B **FUNDRAISING** **PROGRAM**

Electric Rechargeable Systems^{TM*}



TRANSFORMATIONAL GUM HEALTH SYSTEM

- Oral-B® iO™ Series 8 Electric Toothbrush
- PRO-HEALTH™ Gum Detoxify™ Deep Clean Toothpaste
- PRO-HEALTH™ Multi-Protection Clean Mint Mouthwash
- PRO-HEALTH™ Deep Clean Floss

Brush head refills include:



Ultimate Clean (3)



Gentle Care



iO™ GINGIVITIS SYSTEM

- Oral-B® iO™ Series 6 Electric Toothbrush
- PRO-HEALTH™ Gum Detoxify™ Deep Clean Toothpaste
- PRO-HEALTH™ Multi-Protection Clean Mint Mouthwash
- PRO-HEALTH™ Deep Clean Floss

Brush head refills include:



Ultimate Clean



iO™ SERIES 2 DAILY CLEAN SYSTEM

- Oral-B® iO™ Series 2 Electric Toothbrush
- Crest® PRO-HEALTH™ Clean Mint Toothpaste
- Crest® PRO-HEALTH™ Advanced Multi-Protection Mouthwash
- PRO-HEALTH™ Deep Clean Floss

Brush head refills include:



Gentle Care

All fundraising orders must be placed via their Crest® + Oral-B® representative (TAE, SSM, PAR) by emailing CrestOralBSchools.IM@PG.com

*Offer subject to change at any time; while supplies last.

continuing the care that starts in your chair



RECOMMEND PRODUCTS FOR YOUR PATIENTS, FAMILY AND FRIENDS THROUGH THE

Crest + Oral-B FUNDRAISING PROGRAM

Crest® Whitening



**CREST® 3DWHITE DAILY
WHITENING SERUM**
80749555



**CREST® 3DWHITESTRIPS™
PROFESSIONAL EFFECTS**
80702652



**CREST® 3DWHITESTRIPS™
SUPREME PROFESSIONAL
EXCLUSIVE + LIGHT**
80725895

EXPECTED RESULTS

Visibly whiter teeth in 2 weeks

Start to see a whiter smile in 3 days with full results in 20 days.

1 hour for visible results.
Professional-level teeth whitening in just 7 treatments.

MODE OF ACTION & DELIVERY SYSTEM

Highly active micro-droplets with over 20% hydrogen peroxide suspended in a water-resistant base whiten while keeping teeth hydrated.
3% peroxide overall.

10% hydrogen peroxide

11% hydrogen peroxide.

APPLICATION TIME

30 minutes (2x per day)

45 minutes.
Maximum once a day.

1 hour + LED light time.
Maximum once a day.

- All Crest at-home whitening products are enamel safe.
- Results depend on different factors, including a patient's age, how yellow their teeth are, and how often they apply the product.
Most patients see teeth whitening results after just a few uses.
- Duration of results varies, depending on the patient's use of teeth-discoloring beverages like wine, coffee, tobacco, soda, etc.
- Remember, at-home whitening products won't work on caps, crowns, veneers, partials, fillings, laminates, or other similar dental work.

All fundraising orders must be placed via your Crest® + Oral-B® representative (TAE, SSM, PAR) by emailing CrestOralBSchools.IM@PG.com

© 2025 P&G ORAL-23202

continuing the care that starts in your chair

