

the Biodex Gait Trainer 2, and energy expenditure was evaluated using the 6-minute walk test.

Results

Functions and Activities						Participation			Environment
Level walking	Stairs	Ramps, Hills	Uneven ground, Obstacles	Cognitive demand	Metabolic Energy Consumption	Safety	Activity, Mobility, ADLs	Preference, Satisfaction, QoL	Health Economics

Category	Outcomes	Results**		Sig.*
		3R15 (mechanic)	3R60 (hydraulic)	
Level Walking	Velocity and gait cycle increased, while step inequality decreased with the use of 3R60. Weight bearing on the amputated side and ambulation index increased with 3R60			
	Velocity [m/s]	0.7 ± 0.2	0.8 ± 0.02	++
	Gait cycle [s]	0.6 ± 0.02	0.7 ± 0.02	++
	Step Inequality [cm]	4.5 ± 3.0	2.6 ± 1.7	--
	Weight bearing on amputated side [%]	45.3 ± 1.9	47.2 ± 3.4	++
	Ambulation index score [0-100]	72.5 ± 3.2	82.3 ± 4.9	++
Metabolic Energy Consumption	Physiological Cost Index was lower with the 3R60.			
	Physiological Cost Index [beats/min]	34.9 ± 8.8	27.8 ± 7.6	--
Activity, Mobility, Activities of Daily Living (ADLs)	Ambulation activities were completed faster when using 3R60.			
	Sit down on a chair [s]	36.4 ± 8.0	31.4 ± 2.2	-
	Stand up from a chair [s]	32.8 ± 8.1	28.1 ± 5.7	--
	Descending stairs [s]	15.5 ± 1.7	12.1 ± 1.2	--
	Ascending stairs [s]	17.3 ± 1.3	14.3 ± 0.9	--
	Descending slopes [s]	19.6 ± 1.3	12.9 ± 1.2	--
	Ascending slopes [s]	21.0 ± 1.2	13.3 ± 0.9	--
	Walking over obstacles [s]	2.3 ± 0.7	1.1 ± 0.3	--
Safety	Patients could stay balanced much longer and had higher scores with the 3R60.			
	Berg Balance scale [0-56]	21.0 ± 7.4	49.0 ± 9.9	++
Preference, Satisfaction, Quality of Life (QoL)	Subjects felt much more comfortable with 3R60 regarding comfort, cosmetic appearance and functionality.			
	Patient Satisfaction (SATPRO) [0-45]	24.4 ± 4.9	38.7 ± 4.4	++
	The quality of life did not show any differences between the two prostheses, except for less pain with the 3R60. The total score of the Nottingham Health Profile was lower for 3R60.			
	Nottingham Health Profile Energy Level [0-100]	19.8 ± 8.4	8.7 ± 4.6	0
	Pain [0-100]	14.2 ± 4.9	5.1 ± 3.0	--

Category	Outcomes	Results**		Sig.*
		3R15 (mechanic)	3R60 (hydraulic)	
	Emotional Reaction [0-100]	10.2 ± 5.6	8.2 ± 3.6	0
	Sleep [0-100]	10.6 ± 5.6	0.0 ± 0.0	0
	Social Isolation [0-100]	18.4 ± 8.4	17.8 ± 6.6	0
	Physical Abilities [0-100]	17.9 ± 6.1	12.1 ± 6.0	0
	Total Score [0-600]	91.1 ± 28.6	51.9 ± 12.5	--

* no difference (0), positive trend (+), negative trend (-), significant (++/--), not applicable (n.a.)
significance set at $p < 0.05$; trends set at $0.1 > p > 0.05$;

** all Values are given as Mean ± Standard deviation; KJ: Knee Joint;
bold values indicate significant improvement;

Author's Conclusion	"In conclusion, we assumed that our results were parallel with the relevant literature and may stimulate interest in conducting further work in the field of prosthetic training. We thought that our study which evidences that proper selection of prosthetic components according to the patients functional level has positive effects on quality of life, gait pattern and energy expenditure, may outline a framework for a further guideline designed to describe prosthetic prescription for both private and governmental insurance systems in Turkey." (Ülger et al., 2009)
----------------------------	---

Author's Affiliation	¹ Hacettepe University, Faculty of Health Sciences, Department of Physical Therapy and Rehabilitation
-----------------------------	--

© 2025, Otto Bock HealthCare Products GmbH ("Otto Bock"), All Rights Reserved. This article contains copyrighted material. Wherever possible we give full recognition to the authors. We believe this constitutes a 'fair use' of any such copyrighted material according to Title 17 U.S.C. Section 107 of US Copyright Law. If you wish to use copyrighted material from this site for purposes of your own that go beyond 'fair use', you must obtain permission from the copyright owner. All trademarks, copyrights, or other intellectual property used or referenced herein are the property of their respective owners. The information presented here is in summary form only and intended to provide broad knowledge of products offered. You should consult your physician before purchasing any product(s). Otto Bock disclaims any liability related from medical decisions made based on this article summary.