BruxZir® Esthetic Milling Blanks

Instructions for Use

ENGLISH

Manufactured by

Prismatik Dentalcraft, Inc.

(A wholly owned subsidiary of Glidewell Laboratories)

2144 Michelson Drive

Irvine, CA 92612, USA

Within the U.S.: 888-303-3975

Within Canada: 888-278-0414

1. Product Description

BruxZir® Esthetic blanks are zirconia milling blanks used for the production of highly esthetic zirconia dental restorations. The high esthetics and high strength of BruxZir Esthetic restorations make them viable for use in all regions of the mouth. The manufactured dental restorations are made utilizing a CAD/CAM system for design and manufacture. The designed and manufactured dental restorations are then sintered at a high temperature. BruxZir Esthetic blanks are available in white and pre-shaded varieties that can be easily adjusted to match the final shade.

2. Indications for Use

The device is indicated for use by dental technicians in the construction of custom-made all-ceramic restorations for anterior and posterior locations.

3. Contraindications

None.

4. Technical Data Table

As measured according to ISO 6872:2015

Table 1	
Classification Type 2, Class 5	
Average Flexural Strength* (3-Point Bending)	870 MPa

**BruxZir Esthetic exceeds the ISO 6872 requirements of > 800 MPa when the correct sintering profile is used.*

5. Design of Restoration

Failure to observe the restoration wall and connector thickness recommendations as shown in Table 2 may result in fracture of the final restoration. To avoid warping of multiunit and/or large restorations during sintering, milling a sprue-supporting structure may be advisable. For specific instructions on design of the full-contour dental restorations please refer to the applicable CAD manual designated for your CAD system.

Table 2		
Minimum Design Guidelines for Sintered Restorations		
	Minimum Wall Thickness	Minimum Connector Area
Single Units	0.7 mm	N/A
Anterior Bridges	0.7 mm	> 12 mm ²
Posterior Bridges	0.7 mm	> 14 mm ²

6. Milling

After the restorations have been digitally designed, program your milling machine and enter the enlargement factor (EF) or shrinkage for the blank when prompted to do so. The mill operator should refer to the operating instructions for the milling system for the design/programming method of the restoration, for input of the EF, and for milling the zirconia blanks. The specific EF for the blank used can be located on the label affixed to the side of the BruxZir Esthetic blank. Once milled, carefully remove the restorations from the blank.

⚠ Note: Due to its high density, using a common or generic enlargement factor will

result in incorrect fit. Always ensure the EF listed on the BruxZir Esthetic blank is used.

⚠ Note: If the CAD/CAM mill or tool(s) have been used for materials other than zirconia, they must be cleaned thoroughly before using BruxZir Esthetic. Failure to do so may contaminate the restoration, resulting in decreased strength and/or esthetic performance. Refer to the user manual of the milling system and/or tool(s) for further details.

7. Post Milling


Remove all residual dust from the restorations. Ultrasonic cleaning of restorations is recommended but not required. Ensure that only clean distilled or deionized water is used and the restorations are not left in the ultrasonic unit or left in the water for more than 30 seconds. Failure to follow these recommendations may cause unexpected shade or poor esthetics. Completely dry the restorations before sintering. Drying can be accelerated by the use of microwave convection or a heat lamp.

⚠ Note: BruxZir Esthetic restorations do not need to be dipped in any coloring liquids or undergo any shade alteration prior to sintering.

8. Sintering

All BruxZir Esthetic dental restorations must go through a final sintering cycle prior to finishing. Ensure that all sintering accessories are clean and contamination-free. Place the restorations, occlusal side down, in a pure alumina (≥99.9%) sintering tray containing the appropriate-sized (2-mm diameter recommended) zirconia beads (not provided by Prismatik Dentalcraft, Inc.). Ensure that the tray is covered with an alumina lid. The sintering cycle must follow the profile shown below. Failure to do so may result in decreased strength and esthetic performance.

Ensure that the oven is calibrated, in good working condition, and free of any contamination sources. Program your sintering oven according to the Sintering Temperature Chart. Refer to the user manual of your sintering oven for programming information. Allow the oven to cool sufficiently in order to safely remove the restorations. After sintering, carefully inspect the restorations to verify that there are no cracks or chips.

Temperature (°C)	
C1	25
C2	1,100
C3	1,100
C4	1,200
C5	1,450
C6	1,200
C7	1,200
C8	1,475
C9	1,475
C10	1,550
C11	1,550
C12	155

Time (minutes)	
T1	43
T2	60
T3	50
T4	25
T5	1
T6	90
T7	18
T8	5
T9	8
T10	10
T11	90
T12	Off

Sintering Cycle Requirements

- 1st heating rate – Any rate up to 25 C/min to 1100 C
- 1st holding time at 1100 C – 60 minutes
- 2nd heating rate – 2 C/minute to 1200 C
- 3rd heating rate – 10 C/minute to 1450 C
- 2nd holding time at 1450 C – 1 minute
- 1st cooling rate – max cooling down to 1200 C
- 3rd holding time at 1200 C – 90 minutes
- 4th heating rate – 15 C/minute to 1475 C
- 4th holding time at 1475 C – 5 minutes
- 5th heating rate – 10 C/minute to 1550 C
- 5th (Final) holding time at 1550 C – 10 minutes
- Maximum cooldown rate – 15 C/min

⚠ Note: Failure to use this sintering cycle may cause the strength to be lower than 800 MPa.

9. Adjusting the Sintered Restorations

After sintering, if any adjustments to the restorations are necessary, use appropriate low-speed or water-cooled grinding instruments. Ensure that the minimum wall thickness and connector dimensions are maintained (Table 2). To avoid chipping or fracturing, minimize adjustments on thin margins and embrasures.

10. Achieving Final Shade

The BruxZir Esthetic blanks are individually designed to match the indicated VITA® shade after staining and glazing. Because the blanks are monochromatic, some staining and/or shade alterations may be used to match specific prescription requirements and to achieve a gingival to incisal blend. It is recommended to use low-fusing glaze on the restorations (refer to manufacturer's instructions for application guidelines). Using a porcelain or laboratory oven (not provided by PrismaTik), heat the stained and glazed restoration according to the instructions from the manufacturer. Check the quality of the shade, and re-stain if needed. Failure to follow these recommendations may result in an unintended shade.

⚠ Warning: Do not remove the restorations from the oven until they have cooled sufficiently for safe removal.

⚠ Note: The restorations will appear to be darker and/or higher in chroma out of the sintering oven. The restorations are designed to match the VITA shade after a stain/glaze heat cycle.

⚠ Note: Please refer to the operating instructions for the particular oven used for proper use, placement of the product and calibration. Ensure that the sintering oven used is calibrated according to the particular manufacturer's specifications. Contact your sintering oven manufacturer for further instructions.

11. Cementation

No specific cements have been evaluated for BruxZir Esthetic restorations. Industry-standard materials used to cement zirconia restorations may be used. Refer to the user manuals of those specific products for the zirconia cementation protocols.

12. Storage

Before storing BruxZir Esthetic blanks, ensure there is no damage from shipping. Shipping damage may include cracks and/or chips of the blank. The blanks should be stored at room temperature and away from dust and potential contamination sources. When not in use, store the BruxZir Esthetic blanks in the original packaging.

13. Warnings

Do not inhale zirconia dust. Do not get zirconia dust into your eyes or allow it to have prolonged contact with your skin. Personal protective equipment (PPE) such as gloves, masks, and lab coats are recommended in order to avoid possible contact and irritation to skin or eyes during processing.

⚠ Note: To request the SDS for this product, call Glidewell Direct at 888-303-3975.

14. For More Information

For more information, product features, tips and techniques, refer to the technical and training documentation available at bruxzir.com.

15. Questions/Comments

For questions and/or comments or for further instructions regarding the use of the product, please contact Glidewell Direct in the U.S. at 18651 Von Karman Ave., Irvine, CA 92612 or by calling 888-303-3975.

16. Symbol Glossary

Symbol	Reference No. 21 CFR Part 801, Section	Title - Symbol Description
	Sect. 801.109(b)(1)	By Prescription Only – Caution: Federal law restricts this device to sale by or on the order of a dentist or physician.
Symbol	Reference No. EN ISO 15233-1, Section	Title - Symbol Description
	5.1.1	Manufacturer – Indicates the medical device manufacturer.
	5.1.3	Date of Manufacture – Indicates the date (YYYY-MM-DD) when the medical device was manufactured.
	5.1.5	Lot Number – Indicates the manufacturer's batch code so that the lot can be identified.
	5.1.6	Catalogue Number – Indicates the manufacturer's catalogue number so that the medical device can be identified.
	5.1.7	Serial Number – Indicates the manufacturer's serial number so that a specific medical device can be identified.
	5.3.4	Keep Dry – Indicates a medical device that needs to be protected from moisture.
	5.4.3	Consult Instructions For Use – Indicates the need for the user to consult the instructions for use.
	5.4.4	Caution – Indicates the user needs to consult the instructions for use for important cautionary information such as warnings and precautions.

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