



SCAN FOR TIPS, TRICKS, AND TO STOCK UP ON ADDITIONAL WHITENING PENS!

TRYBYTE.COM.AU/PAGES/BRIGHTBYTE-PRO

SHOW OFF THOSE PEARLY WHITES (AND YOUR TEETH SPA ROUTINE) BY TAGGING

@BYTEOFFICIAL\_AU



## BrightByte™ PRO

WIRELESS WHITENING KIT WITH MULTI-SPECTRUM  
LED ACCELERATOR LIGHTS

REF BBPRO2



Manufactured by  
Nanchang Dental Bright Technology Co., Ltd  
Building 182 Jiahai Industrial Park, High-Tech Development Zone,  
Nanchang, Jiangxi Province, China, 330096

Distributed by  
Byte Au Pty Ltd, 8/226 Condamine Street, Manly Vale NSW 2093

Made in China

Questions? Call 1800 298 363 or go to [patient.byte.com/support](http://patient.byte.com/support)

PRT.BBPRO2-001.0

### Ready to give your oral wellness routine the ultimate spa treatment?

The BrightByte™ Pro is an on-the-go LED whitening kit for achieving your brightest smile. Our exclusive mouth tray uses cutting-edge, low-temperature LED blue light technology, paired with our whitening gel, to remove stains efficiently—or combine your treatment with a gentle purple LED light for a unique whitening experience. Recover after whitening by using the gentle massage feature, or try the soothing red LED light which can help to calm and rejuvenate your teeth and gums.

Use the whitening treatment for 16 minutes a day, and you'll quickly start to see a brighter smile. Perfect for regular whitening maintenance and customising your oral wellness habits.

## HOW TO USE

### STEP 1

Rinse the BrightByte Pro mouth tray and your mouth before use. Dry your teeth to prevent saliva from diluting the gel adherence. Look at yourself in the mirror and compare your teeth to the shade guide (included), then record what stage you're on. Twist open the whitening pen until the gel flows onto the brush. Apply approx. 0.5 ml to 0.7 ml (about a pea-sized amount) evenly to all teeth.

### STEP 2

Put the mouth tray into your mouth and gently bite down with the power button facing up.

### STEP 3

Set a timer to whiten for up to 16 minutes total. When the time is up, remove the mouth tray and rinse the tray and your mouth with water. Compare your teeth with the shade guide again and record your results. Do not eat or drink anything (except for water) for 30 minutes after completing your whitening treatment.

### STEP 4 (OPTIONAL)

If you're feeling any sensitivity on your teeth, use Mode 3 for a blue light and gentle massage or Mode 4 for a red light which may help soothe teeth after treatment. 5 minutes total per day.

### QUICKSTART GUIDE:

Mode	Function	Recommended usage
1. Purple light	Combination whitening	16 min max per day with whitening gel
2. Blue light	Basic whitening	
3. Blue light + massage	Massage treatment	5 min max per day
4. Red light	Recovery treatment	

## SAFETY DIRECTIONS

Avoid contact with eyes, skin and wash hands after use.

### FIRST AID

If in eyes, hold eyelids apart and flush the eye continuously with running water. Continue flushing until advised to stop by a Poisons Information Centre. If swallowed, do NOT induce vomiting. If skin or hair contact occurs, remove contaminated clothing and flush skin and hair with running water. For advice, contact a Poisons Information Centre (e.g. phone Australia 13 11 26; New Zealand 0800 764 766) or a doctor (at once).

## HOW TO TOGGLE BETWEEN MODES

- Press the power button once to turn on the device in Mode 1. Red and blue lights will illuminate to create a purple hue.
- Press the power button a second time to put the device in Mode 2. Blue light will illuminate.
- Press the power button a third time to put the device in Mode 3. Blue light will illuminate and a gentle massage will begin.
- Press the power button a fourth time to put the device in Mode 4. Red light will illuminate.
- Press the power button a fifth time to turn off the device.

## CAUTION Irritant.

- Do not apply the gel directly on your gums or on your inside lips, as this may create discomfort and irritation in the applied area. If you experience discomfort, rinse with water.
- Not suitable for veneers, ceramic teeth, or other dental restorations.
- Not suitable for tetracycline stained teeth, decayed teeth, defective enamel, exposed dentin, damaged teeth, or teeth discoloration caused by lesion or medication.
- Not suitable for children or anyone that is pregnant (or nursing).
- Do not use gel pens after the expiration date.
- Replace cap after use. Avoid sunlight; store in a cool, dry place.
- DO NOT USE BrightByte Pro gentle massage feature in place of HyperByte® during aligner treatment.
- Do not use this light if you are being treated with photodynamic therapy by your doctor, are taking light sensitivity medicines like aminolaevulinic or porphyrin containing medicines or have a sensitivity to light.

## INGREDIENTS

### Whitening Gel

Carbamide (Urea) Peroxide, Glycerin, Propylene Glycol, Aqua, Carbomer, Cellulose Gum, PVP, Menthol, Sodium Hydroxide.

## FAQ

- Q:** My BrightByte Pro won't turn on?  
**A:** Charge the device for more than 2.5 hours, press the power button, and check to see if the LED light is functioning normally.
- Q:** What types of stains does BrightByte Pro work on?  
**A:** Using the whitening gel for up to 16 minutes daily will help to improve the appearance of surface level stains for a white smile. Treatment is enhanced with the wireless LED accelerator light.
- Q:** Is the blue light technology safe?  
**A:** The device's LED light wavelengths do not cause any irritation or damage to tooth enamel or gingiva.
- Q:** How long will the whitening effect last?  
**A:** Daily habits and diet can affect tooth appearance over time, and whitening maintenance may be required to keep your smile white.