

Fixodent

4 EASY STEPS



CLEAN & DRY

1

Clean and dry your denture before you apply the adhesive cream.



APPLY

2

Apply adhesive in thin strips as shown or less. You might need a few tries to find out the ideal amount for your individual denture.



Upper



Lower



Partial

Maximum size of each strip

4x16.5mm



INSERT & HOLD

3

Insert denture and hold it briefly in place. If product oozes off denture in mouth, you are using too much.



RINSE & REMOVE

4

To remove your denture, rinse your mouth with warm water and rock denture gently until loose.

UK IE ZA ADHESIVE CREAM for full and partial dentures.
IMPORTANT:

- ✗ **DO NOT** use more product than shown.
- ✗ **DO NOT** use product more than once a day. With proper use, this tube should last 7 - 8 weeks for 70g and 5 - 6 weeks for 47g.
- ✗ **DO NOT** use excess product for poorly fitting dentures.
- ✗ **DO NOT** use more than directed. Contains zinc, consult a doctor if using other products containing zinc.

Doesn't feel right? Please consult your dentist regularly to ensure you have properly fitting dentures. Poorly fitting dentures may impair your health. Once opened, use within 6 months.

WARNING: Keep out of reach of children. Do not use in the case of hypersensitivity to ingredients. In the unlikely event of a serious suspected problem with this product consult your dentist or doctor and report this to the manufacturer and national health authority (UK: MHRA; IE: HPRA). Keep nozzle and cap dry to avoid clogging. If tip becomes clogged, use a pin to reopen.

COMPLETE No Flavour INGREDIENTS: Calcium/Zinc PVM/MA Copolymer (33%), Paraffinum Liquidum, Petrolatum, Cellulose Gum (20%), Silica, BHT, Tocopherol.

COMPLETE Fresh Flavour INGREDIENTS: Calcium/Zinc PVM/MA Copolymer (33%), Paraffinum Liquidum, Petrolatum, Cellulose Gum (20%), Silica, Mentha Piperita Oil, CI 45410, Menthol, Limonene, Beta-Caryophyllene, Pinene, Aroma, BHT, Tocopherol.

Made in Germany for

P&G Tech. Centres Ltd, KT13 OXP, UK

UK Procter & Gamble UK, Weybridge, Surrey KT13 OXP. Tel: 0800 731 3086.

IE 1800 409 236 **ZA** Distributed by P&G South African Trading (Pty) Ltd, 15 Alice Lane, Sandton, 2196 Gauteng, South Africa. Consumer Care Line 0860 112 188 (Sharecall charged at local rates). www.pg.com/contactus.

Issue: 2017-11

CE 0197



www.pg.com