Gleo Fluoride Varnish 5% Sodium Fluoride Varnish

Product Description
Gleo fluoride varnish is a flavored, xylitol-sweetened, 5% sodium fluoride in a resin carrier. For all products described, carefully read and understand all instructions and SDS information prior to use.

Delivery Form(s)

- 10 ml tube
 - Application brush

Indications for Use/Intended Purpose

Gleo fluoride varnish is 5% sodium fluoride in a varnish carrier which produces a mechanical occlusion of the dental tubules in the treatment of tooth hypersensitivity.

Contraindications

- Do not use on patients with ulcerative gingivitis or stomatitis.
- Do not use on patients with a known sensitivity to ingredients including colophonium (kolophonium) and/rosin. If allergic reaction, dermatitis, or rash occurs, discontinue use.
- For patients or users with allergy concerns, refer to product allergen document available at ultradeut.com. If allergic reaction is observed, rinse exposed and thoroughly with water and have the patient consult their physician.

Warnings and Precautions

- Patient, clinician, and assistant must wear protective eyewear.
- In rare instances, edematous swelling has been reported after extensive fluoride varnish treatment and dyspnea has been reported in asthmatic children. If this occurs, thoroughly brush teeth to remove the varnish film and have patient contact their medical provider.
- Some patients may experience nausea after application. If nausea persists, thoroughly brush to remove the varnish film.
- Do not use with other topical fluoride treatments the same day and instruct patient to refrain from using oral fluoride supplements for four days after treatment.
- Only use clean and sanitized instrument to coat upon the tube.

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