Gishy Goo™
Silicone Bracket Relief Aid
EN

Product Description
Gishy Goo™ Silicone Bracket Relief Aid is a flavored, thixotropic, silicone elastomer that holds firmly to brackets and wires. It comes in a syringe with two barrels. One barrel has a white activator material. The other has a light yellow base material. A pea-sized amount of Gishy Goo material is dispensed by the patient into their hand, then mixed and applied to bracket or wires at the point of irritation. The product sets in two minutes and stays in place for up to 12 hours.

Gishy Goo material can be used on patients of all ages and conditions requiring a bracket relief aid other than those with known sensitivities to silicone. Any restrictions regarding patient populations can be found below in the Contraindications section.

For all products described, carefully read and understand all instructions and SDS information prior to use.

Delivery Form(s)
- 5mL dual barrel syringe

Indications for Use/Intended Purpose
Gishy Goo Silicone Bracket Relief Aid provides mechanical relief and a protective barrier to prevent irritation of soft tissue caused by orthodontic appliances. Gishy Goo Silicone Bracket Relief Aid is applied to points of irritation on brackets, wires and fixed/functional appliances.

Contraindications
- Do not use on patients with known allergies to silicones.
- For patients or users with allergy concerns, refer to product allergen document available at ultradent.com/resources/product-allergens. If allergic reaction is observed, rinse exposed area thoroughly with water and have the patient consult their physician.

Warnings and Precautions
- Children under the age of 12 should be supervised or assisted by parents when using and applying Gishy Goo material. Gishy Goo material may represent a choking hazard to children 4 years of age and younger.
- Take Gishy Goo material off brackets before sleeping, eating, or drinking to avoid the chance of choking.
- Do not eat; large amounts may present a choking hazard.
- Do not handle using latex gloves.
- The maximum amount of Gishy Goo material for each sore spot is about the size of a pea. If you squeeze out too much, throw it away and try again.
- You should be mixing equal amounts of white and colored Gishy Goo material. If you don't have about the same amounts, throw away and try again.
- Do not store Gishy Goo material near heaters or leave it in hot cars, or it might not set properly.
- Do not squeeze or stick Gishy Goo material on clothes or fabric, as it will not be removable.

Stepwise Instructions

Preparation
1) Remove cap by rotating counterclockwise.

Use
1) Squeeze out a pea-sized amount onto finger or hard surface. Be sure to get the same amount of both white and colored material. If only one color comes out, throw that material away and try again.
2) Mix by hand for about 10 seconds.
3) Make sure area is fairly dry and place mixed Gishy Goo material on the brackets at the point of irritation. Gishy Goo material takes about 2 minutes to set.
4) After use, put the cap back on. Make sure cap is on tightly.

Cleanup
1) To remove, use a toothpick or a toothbrush to push it off.
2) Gishy Goo material may remain in place for up to 12 hours. Remove all Gishy Goo material from brackets before going to sleep, eating, or drinking.

Storage and Disposal
- Store at room temperature.
- Dispose of waste according to local rules, guidelines, and regulations.

Miscellaneous Information
Report any serious incident to the manufacturer and the competent authority.
ultradent.com / 800.552.5512 / 801.572.4200

Consult healthcare professional for any product related questions.

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