

Point of Care Testing Accuchek Inform Post Test

NAME _____ EMPLOYEE # _____ DATE _____ UNIT _____

Directions: Circle the correct response.

| | | | |
|---|---|--|---|
| T | F | | 1. If the patient's fingers are cyanotic or edematous; the bedside glucose test result will not be valid. Perform the bedside test on a venous or arterial specimen instead; or order a Lab draw glucose. |
| T | F | | 2. Control solutions are stable indefinitely. |
| T | F | | 3. The best results are obtained by wiping away the first drop of blood and using the second drop of blood for testing. |
| T | F | | 4. The patient financial (admit) number must be entered for the test to record to the proper patient account and be billed. |
| T | F | | 5. If the yellow visual window is not completely filled, a second drop may be applied within 15 seconds of the first drop. |
| T | F | | 6. The glucose meter must be docked to recharge and download results to the patient chart. |
| T | F | | 7. Laboratory confirmation is required for ALL glucose tests performed at the bedside. |
| T | F | | 8. It is helpful to know the patient's Hematocrit when interpreting the validity of the Accuchek glucose results. |
| T | F | | 9. The Accuchek glucose meter is cleaned with a wrung out alcohol pad after each patient use. |
| T | F | | 10. Comment codes are required for every test. |
| T | F | | 11. After removing a test strip from the closed vial, it must be used within 10 minutes. |

Pass (100%) _____ Fail _____

Certified Trainer Signature: _____ Date: _____