Which CGM will you choose?

Continuous Glucose Monitors Comparison

This comparison tool will help answer common questions to help make choosing one of the two most common brands of CGMs easier. **Contact your insurance provider to determine if a CGM is covered under your insurance plan and if your insurance has a preferred brand.**

Features	FreeStyle Libre 2/3 from Abbott	G6 or G7 from Dexcom
How long does a sensor last?	14 days, no grace period	G6: 10 day, no grace period G7: 10 day, +12 hours beyond expiration (grace period)
Is it compatible with iPhone and Android devices?	YES	YES
I don't have a smartphone, is there a reader I can use with the CGM?*	YES	YES
Is a separate transmitter required? If yes, how long does the transmitter last?	NO	G6: YES; 90 days G7: NO
What is the communication type (i.e. NFC, Bluetooth)?	Libre 2: NFC to start and scan the sensor, Bluetooth for alarms. Libre 3: NFC to start the sensor, Bluetooth for streaming glucose data and alarms	Bluetooth

*If you use a standalone receiver, you will need to upload your receiver to Libreview (FreeStyle Libre) or Clarity (Dexcom) before your doctor appointment, or bring it with you to your appointment if you do not have computer and internet access at home.

Applying the CGM	FreeStyle Libre 2/3 from Abbott	G6 or G7 from Dexcom
What is the size of the sensor?	Libre 2: 1.38in. wide x 0.20in. tall Libre 3: 0.83in. wide x 0.11in. tall	G6: 1.2in. wide x 1.8in. long G7: 1.07in. wide x 0.94in. long
Is the applicator pre-attached to sensor?	Libre 2: NO Libre 3: YES	G6 or G7 : YES
Where on the body can the sensor be worn?	Back of upper arm	G6: Abdomen G7: Back of upper arm
What is the warm-up period (the period of time between applying the CGM and being able to read your glucose level)?	1 hour	G6: 2 hour G7: 30 min

Using the CGM	FreeStyle Libre 2/3 of from Abbott	G6 or G7 from Dexcom
What required alarms does the CGM have?	Urgent low alarm (<55mg/dL)	Urgent low alarm (<55mg/dL)
What optional alarms does the CGM have?	High glucose, low glucose, and signal loss alarms	High & low glucose alerts, predictive urgent low soon
Is calibration required?	NO	NO. Optional calibration is available.
Can I share my CGM data with my care team?	Via LibreView website	Via Clarity app and website



Support	FreeStyle Libre 2/3 from Abbott	G6 or G7 from Dexcom
What is the replacement policy for faulty sensors or sensor that fall off prematurely? How can I obtain a replacement sensor?	Answer the guided questions and complete the request form: <u>michmed.org/RnRMW</u>	Complete a Dexcom Product Support Request at <u>michmed.org/8zJ4G</u>
Where can I go to find educational materials about this CGM? MCT2D CGM Resources michmed.org/JYQZd	MyFreeStyle Free guided support for new FreeStyle users FreeStyle Libre 2 User Guide: michmed.org/38xRx FreeStyle Libre 3 User Guide: michmed.org/wr8Xr	Dexcom Care Virtual Program – <u>michmed.org/5k8WW</u> register to meet virtually with a Certified Diabetes Care and Education Specialist (CDCES) Dexcom Training, User Guides, & Support Materials: <u>michmed.org/</u> <u>W44jZ</u>
What live (phone or text) customer support is offered?	Customer Care: 1-855-632-8658 Also available via Facebook Messenger.	Dexcom Customer Support: 1-888-738-3646 Online Chat (8am-8pm EST): <u>michmed.org/kqqgK</u>

Additional CGM Resources



Free Trial CGM



FreeStyle Libre 3 or 2 CGM michmed.org/XYwD7







Dexcom G7 CGM

michmed.org/xANnx



Visit GoodRx.com or download the app to find coupons for CGMs