



ottobock.

connectgrip.

Real-time visualization with
myoelectric signal feedback.

The **connectgrip** app provides a user-friendly interface for users and professionals to control and personalize myoelectric arm and hand prostheses. It enables easy adjustments, real-time signal visualization, and individual motion recording for movements, joint switching and unintended movements. Users achieve improved control, while professionals benefit from streamlined prosthetic setup features, improving both functionality and comfort.




Download on the
App Store



GET IT ON
Google Play

connectgrip.

Optimal flexibility in customizing users' myosignals.



- Personalization with advanced signal processing and independent adjustments for everyday use.
- User friendly interface for easy daily adjustments.
- Real-time visualization with immediate myoelectric signal feedback.
- Individual motion recording to minimize unintended movements.
- Prosthetic setup and documentation streamlining the fitting process for healthcare professionals.

Get started in three easy steps

1

Install: Download the **connectgrip** app from Google Play or the Apple Store.

2

Log in: Open the **connectgrip** app and log in with your MyOttobock Pro account via the user icon or in the tab bar under 'Menu'. If you do not yet have a MyOttobock account, please register via the 'Register here' link on the login screen. The **connectgrip** app automatically directs you to the registration of a 'MyOttobock Pro' account. This is the version of MyOttobock for prosthetists and orthotists.

3

Activate: If you have already completed a product certification via MyLearning in advance, this will be automatically added to your MyOttobock Pro account, and you can use the **connectgrip** app in specialist mode to set up the product. Alternatively, you can request an activation for the product via the user icon or the 'Menu / Certification' tab bar. To do this, follow the instructions on the screen.

Manufactured by:
Otto Bock Healthcare Products GmbH
Brehmstraße 16 Wien, Wien, 1110 Austria

Distributed by:
Otto Bock HealthCare LP
3820 Great Lakes Drive
Salt Lake City, UT 84120

<https://www.ottobock.com/en-us/professionals/prosthetics/connectgrip>

US: P 800 328 4058 | F 800 962 2549
CA: P 800 665 3327 | F 800 463 3659