

ottobock.



This single RevoFit Lamination Kit with the Click Reel allows you to fabricate user-controlled adjustability anywhere on a device.

The Click Reel innovation creates adjustable devices with a lower streamlined profile, micro/macro adjustment that incrementally tightens and loosens with increased durability.



Information material

646D1752=EN_INT	Product Brief RevoFit
-----------------	-----------------------

Product example



21Y164=1 **RevoFit®**

The RevoFit® closure and adjustment mechanism with Click Reel® technology offers a better fit and optimized comfort of the prosthesis. Thanks to the Click Reel® system, micro-adjustment of the socket is possible.

Kit includes:

1x Click Reel, 1x Lamination Collar, 2m HD Lace, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 1.8m Lamination Tubing

User Weight Limit: 300lb (135kg)

Technical data

Article number	Quantity unit
21Y164=1	Set

21Y164=10 **RevoFit®**

The RevoFit Lamination Fabricators 10-Pack with the Click Reel is recommended for central fabricators and practices that create large volumes of adjustable devices.

Technical data

Article number	Quantity unit
21Y164=10	Set



21Z18=1 **Lace Kit for RevoFit**

2 m of HD (High Durability) Lace and Click Reel components for RevoFit Lamination Kits with the Click Reel.

Technical data

Article number
21Z18=1





21Z19=1 Click Reel, Lace Kit for RevoFit

Replacement Click Reel and 2 m of HD (High Durability) Lace for RevoFit Kits only. This is not for replacing a BOA® dial. Please see 21Z37=1 for replacing a BOA® dial.

Technical data

Article number

21Z19=1



21Z15=1 Lamination Collar

Replacement Lamination Collar for RevoFit & RevoLock Lamination Kits with the Click Reel only.

Technical data

Article number

21Z15=1



21Z16=1 Lamination Tube – 1.8 m

Replacement Lamination Tube (1.8 m) for all Revo Kits with the Click Reel and BOA® dials.

Technical data

Article number

21Z16=1



21Z20=1 Spool with Lace for RevoFit

25 m spool of HD (High Durability) Lace assures you will always have replacement lace on hand.

Cost effective, sustainable and less waste.

Compatible with both the RevoFit Lamination Kit with the Click Reel and RevoFit Lamination Kit with BOA® dial.

Technical data

Article number

21Z20=1



21Z22=1 Slider Part Kit

Replacement Slider and Lock Wire for dynamic AFOs, large gap prostheses or RevoLock applications.

Technical data

Article number

21Z22=1



21Z23=1 Slider Lamination Kit

Lamination Kit for fabricating new dynamic AFOs, large gap prostheses or replacement for RevoLock Kits.

Technical data

Article number

21Z23=1



21Z17=1 Grip Ring for Click Reel

Provide a soft-touch feel and 10% increased turning power to the existing Click Reel. Great for patients with dexterity issues.

Technical data

Article number

21Z17=1



21Z17=2 Grip Ring for Click Reel

Provide a soft-touch feel and 20% increased turning power to the existing Click Reel. Great for patients with dexterity and hand strength issues.

Technical data

Article number

21Z17=2



21Z43=1 Click Care Kit™

The Click Care Kit enables you to service your clients in the office or on the road with everything readily accessible for maintaining, repairing or replacing Click components.

Contents:

1x Surface Collar Kit, 2x HD Lace Kits, 2x HC Lace Kits, Flathead Screwdriver, Dental Pick, Spectra Scissors, Needle Nose Pliers & Reel Lubricant. Extra space is allocated for a Click® Reels and a lighter, which can be purchased separately.

Technical data

Article number

21Z43=1

This single RevoSurface Tool Kit with the Click Reel allows you to test adjustable designs and fit of check sockets, 3D printed devices, and on orthoses. It can also be used to solve fit issues on existing rigid devices.

The Click Reel innovation creates adjustable devices with a lower streamlined profile, micro/macro adjustment that incrementally tightens and loosens with increased durability.



Information material

646D1754=EN_INT	Product Brief RevoSurface
-----------------	---------------------------

Product example



21Y167=1 **RevoSurface®**

The RevoSurface® system with Click Reel® technology can be used in a very wide range of applications – wherever contraction and opening of different components is required.

Kit includes:

1x Click Reel, 2x Surface Collar,
2m HC Lace, 1.8m Lamination Tubing, 1x Metal
Lace Feeder, 3x Plastic Lace Feeders, 1x Reel
Tool, 6x Surface Guides, 3x Releasing Guides,
4x Chicago Screw assemblies

User Weight Limit: 300 lb (135 kg)

Technical data

Article number	Quantity unit
21Y167=1	Set



21Z13=1 Indicator Lace Kit

2 m of HC (High Capacity) Lace and Click Reel components for RevoSurface Tool Kit, RevoLock Align. and RevoLock 4-Hole Kits with the Click Reel.

Technical data

Article number	21713=1
-----------------------	---------



21Z12=1 Click® Reel and Lace Kit

Replacement Click Reel and 2 m of HC (High Capacity) Lace for RevoLock Align, RevoLock 4-Hole and RevoSurface Tool Kits only.

This is not for replacing a BOA® dial. Please see 21Z37=1 for replacing a BOA® dial.

Technical data

Article number	21712=1
-----------------------	---------



21Z9=1 Surface Guide 6-Pack

The 6-Pack of Surface Guides allows a practice to have on hand additional guides for easy replacement and configuration changes. For RevoSurface Tool Kits with the Click Reel only.

Technical data

Article number
21Z9=1



21Z14=1 Surface Collar

2x Click Reel mounting collar. For surface mount only. Not for use with lamination.

Technical data

Article number
21Z14=1



21Z10=1 Surface Releasing Tab 3-Pack

3-Pack Surface Releasing Tabs for converting a device from a fixed lace path to a releasing lace path to fully open the device or to replace existing Surface Releasing Tabs.

Technical data

Article number
21Z10=1



21Z16=1 Lamination Tube – 1.8 m

Replacement Lamination Tube (1.8 m) for all Revo Kits with the Click Reel and BOA® dials.

Technical data

Article number
21Z16=1



21Z11=1 Indicator Lace – 25 m Spool

25 m spool of HC (High Capacity) Lace assures you will always have replacement lace on hand.

Cost effective, sustainable and less waste.

Used in the RevoSurface and RevoLock Kits with the Click Reel only.

Technical data

Article number

21Z11=1



21Z17=1 Grip Ring for Click Reel

Provide a soft-touch feel and 10% increased turning power to the existing Click Reel. Great for patients with dexterity issues.

Technical data

Article number

21Z17=1



21Z17=2 Grip Ring for Click Reel

Provide a soft-touch feel and 20% increased turning power to the existing Click Reel. Great for patients with dexterity and hand strength issues.

Technical data

Article number

21Z17=2



21Z43=1 Click Care Kit™

The Click Care Kit enables you to service your clients in the office or on the road with everything readily accessible for maintaining, repairing or replacing Click components.

Contents:

1x Surface Collar Kit, 2x HD Lace Kits, 2x HC Lace Kits, Flathead Screwdriver, Dental Pick, Spectra Scissors, Needle Nose Pliers & Reel Lubricant. Extra space is allocated for a Click® Reels and a lighter, which can be purchased separately.

Technical data

Article number

21Z43=1

This kit offers the push to engage and pull to release the interface of the BOA® Dial in a RevoLock Lanyard solution. It includes a shorter and smaller diameter distal housing which accepts the Quick Connect Snap for easy connection to the liner.



21Y168=1 **RevoLock® Lanyard V3 with BOA® Dial**

This kit offers the push to engage and pull to release interface of the BOA® Dial in a Revo-Lock Lanyard solution. It includes a shorter and smaller diameter distal housing which accepts the Quick Connect Snap for easy connection to the liner.

Kit includes:
1x BOA® Lanyard Dial, 1x Lamination Collar, 1x 0.5 m Lamination Tubing, 2x LE Quick Connect Inserts, 1x Quick Connect Snap, 1x Align Housing + Fabrication Post, 2x 1.22 m HC Lace, 1x Metal Lace Feeder

User Weight Limit: 220 lb (100 kg)

Technical data
Article number
21Y168=1



21Z28=1 **RevoLock® Lanyard Spare parts set (Port A)**

Replacement 1x BOA® Lanyard Dial, 1.22 m of HC Lace (High Capacity), 1x Metal Lace Feeder. For use with Port A Lace exit.

Compatible only with RevoLock Lanyard V2 and V3 Kits.

Technical data
Article number
21Z28=1=1



21Z28=2 **RevoLock® Lanyard Spare parts set (Port C)**

Replacement 1x BOA® Lanyard Dial, 1.22 m of HC Lace (High Capacity), 1x Metal Lace Feeder. For use with Port C Lace exit.

Compatible only with RevoLock Lanyard V2 and V3 Kits.

Technical data
Article number
21Z28=2



21Z29=1 **RevoLock® Lanyard HC Lace Kit**

1.22 m of HC Lace (High Capacity), 1x Metal Lace Feeder.

Compatible only with RevoLock Lanyard V2 and V3 Kits

Technical data
Article number
21Z29=1

**Information material**

646D1753=EN_INT	Product information for technicians/qualified personnel
-----------------	---

Product example

21Y165=1 RevoLock®

This single RevoLock 4-Hole Kit with the Click Reel is for standard alignment fabrication for lower extremity sockets and comes with a standard 25 mm x 65 mm distal housing and integrated 4-Hole plate.

The Click Reel innovation offers increased durability and power, a lower streamlined profile and micro/macro adjustment that incrementally tightens and loosens. The Slider enhances safety for limbs changing in position or fluctuating in volume. RevoLock Tool Form required for fabrication of RevoLock 4-Hole kit.

Closure and adjustment mechanism to secure liner in a prosthetic socket

Kit includes:

1x Click Reel,
1x Lamination Collar, 1x Slider Lamination Kit, 2x RevoLock LE Quick Connect Inserts, 1x RevoLock LE Quick Connect Snap, 2 m HC Lace, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 0.5 m Lamination Tubing, 1x 4-Hole Puck, 1x Tool Form Plug

User Weight Limit: 300 lb (135 kg)

Technical data

Article number	Quantity unit
21Y165=1	Set

**Information material**

646D1753=EN_INT	Product information for technicians/qualified personnel
-----------------	---

Product example

21Y166=1 RevoLock®

This single RevoLock Align Kit with the Click Reel is for transferable alignment fabrication for both upper and lower extremity sockets.

The Click Reel innovation offers a lower streamlined profile, micro/macro adjustment that incrementally tightens and loosens with increased durability. The RevoLock Slider enhances safety for limbs changing in position or fluctuating in volume.

Flexible positioning of the reel housing via individual selection of the course of the pull cord

Kit includes:

1x Click Reel, 1x Lamination Collar,
1x Slider Lamination Kit, 2x RevoLock LE Quick Connect Inserts, 1x RevoLock LE Quick Connect Snap, 2 m HC Lace, 1x Metal Lace Feeder, 3x Plastic Lace Feeders, 1x Reel Tool, 0.5 m Lamination Tubing, 1x Align Housing + Fabrication Post

User Weight Limit: 300 lb (135 kg)

Technical data

Article number	Quantity unit
21Y166=1	Set



21Z4=1 RevoLock® Tool Form

The 25 mm x 65 mm RevoLock Fabrication Tool Form is a required, reusable fabrication tool for RevoLock 4-Hole Kits.

Technical data

Article number
21Z4=1



21Z12=1 Click® Reel and Lace Kit

Replacement Click Reel and 2 m of HC (High Capacity) Lace for RevoLock Align, RevoLock® 4-Hole and RevoSurface Tool Kits only.

This is not for replacing a BOA® dial. Please see 21Z37=1 for replacing a BOA® dial.

Technical data

Article number
21Z12=1



21Z13=1 Indicator Lace Kit

2 m of HC (High Capacity) Lace and Click Reel components for RevoLock Align, RevoLock 4-Hole, and RevoSurface Tool Kits with the Click Reel.

Technical data

Article number
21Z13=1



21Z22=1 Slider Part Kit

Replacement Slider and Lock Wire for dynamic AFOs, large gap prostheses or RevoLock applications.

Technical data

Article number
21Z22=1



21Z23=1 Slider Lamination Kit

Lamination Kit for fabricating new dynamic AFOs, large gap prostheses or replacement for RevoLock Kits.

Technical data

Article number

21Z23=1



21Z15=1 Lamination Collar

Replacement Lamination Collar for RevoLock Align, RevoLock 4-Hole, RevoLock Lanyard V3 & RevoFit Lamination Kits.

Technical data

Article number

21Z15=1



21Z16=1 Lamination Tube – 1.5 m

Replacement Lamination Tube (1.8 m) for all Revo Kits with the Click Reel and BOA® dials.

Technical data

Article number

21Z16=1



21Z5=1 RevoLock® LE Quick Connect Insert

Replacement RevoLock LE Quick Connect Insert for connection to lower extremity liners. Compatible with RevoLock LE Quick Connect Snap (21Z6=1) found in the RevoLock Align, RevoLock 4-Hole and RevoLock Lanyard V3 Kits. Often purchased for outfitting additional patient liners.

Technical data

Article number

21Z5=1



21Z6=1 RevoLock® LE Quick Connect Snap

Replacement RevoLock LE Quick Connect Snap connects the HC (High Capacity) Lace to the liner via the RevoLock LE Quick Connect Insert (RP-302).

Technical data

Article number

21Z6=1



21Z7=1 RevoLock® Tool Form Plug 5-Pack

The RevoLock Tool Form Plug 5-Pack is to ensure fabricators have extra plugs on hand.

Technical data

Article number

21Z7=1



21Z8=1 RevoLock® UE Fixed Insert

RevoLock Fixed Threaded Insert (M6) for upper extremity and pediatric liners. Compatible with RevoLock 4-Hole, RevoLock Align and RevoLock Lanyard V3 Kits. Best for patients who cannot operate the RevoLock Quick Connect Snap or for suction systems needing an airtight seal.

This part is NOT releasable and does not connect to the RevoLock LE Quick Connect Snap.

Technical data

Article number

21Z8=1



21Z21=1 RevoLock® LE Fixed Insert

RevoLock Fixed Threaded Insert (M10) for lower extremity liners only. Compatible with RevoLock 4-Hole, RevoLock Align and RevoLock Lanyard V3 Kits. Best for patients who cannot operate the RevoLock Quick Connect Snap or suction systems needing an airtight seal.

This part is NOT releasable and does not connect to the RevoLock LE Quick Connect Snap.

Technical data

Article number

21Z21=1



21Z11=1 Indicator Lace – 25 m Spool

25 m spool of HC (High Capacity) Lace assures you will always have replacement lace on hand.

Cost effective, sustainable and less waste.

Used in the RevoSurface and RevoLock Kits with the Click Reel only.

Technical data

Article number

21Z11=1



21Z17=1 Grip Ring for Click Reel

Provide a soft-touch feel and 10% increased turning power to the existing Click Reel. Great for patients with dexterity issues.

Technical data

Article number

21Z17=1



21Z17=2 Grip Ring for Click Reel

Provide a soft-touch feel and 20% increased turning power to the existing Click Reel. Great for patients with dexterity and hand strength issues.

Technical data

Article number

21Z17=2



21Z43=1 Click Care Kit™

The Click Care Kit enables you to service your clients in the office or on the road with everything readily accessible for maintaining, repairing or replacing Click components.

Contents:

1x Surface Collar Kit, 2x HD Lace Kits, 2x HC Lace Kits, Flathead Screwdriver, Dental Pick, Spectra Scissors, Needle Nose Pliers & Reel Lubricant. Extra space is allocated for a Click® Reels and a lighter, which can be purchased separately.

Technical data

Article number

21Z43=1



514Z22=XS-1 QuickFit® Lacer

The extra small size QuickFit Lacer is designed to close a small region of an orthosis, a pediatric sized orthosis, or pediatric footwear up to size 35 EU.

Comes pre assembled with a low-profile dial and 4 guides. QuickFit Rivets sold separately, but are required to install the QuickFit Lacer kits. Lacers are sold individually. For a pair of shoes, two kits are required.

Technical data
Article number
514Z22=XS-1



514Z22=XS-10 QuickFit® Lacer

The extra small size QuickFit Lacer is designed to close a small region of an orthosis, a pediatric sized orthosis, or pediatric footwear up to size 35 EU.

Comes pre assembled with 10 extra small QuickFit Lacer kits and 100 QuickFit Rivets. For a pair of shoes, two kits are required.

Technical data
Article number
514Z22=XS-10



514Z22=S-1 QuickFit® Lacer

The small size QuickFit Lacer is designed to close a small region of an orthosis, a small sized orthosis, or footwear size 36 - 41 EU.

Comes pre assembled with a low profile dial and 6 guides. QuickFit Rivets sold separately, but are required to install the QuickFit Lacer kits. Lacers are sold individually. For a pair of shoes, two kits are required.

Technical data
Article number
514Z22=S-1



514Z22=S-10 QuickFit® Lacer

The small size QuickFit Lacer is designed to close a small region of an orthosis, a small sized orthosis, or footwear size 36 - 41 EU.

It comes with 10 Small QuickFit Lacer kits and 100 QuickFit Rivets. For a pair of shoes, two kits are required.

Technical data
Article number
514Z22=S-10



514Z22=M-1 QuickFit® Lacer

The medium size QuickFit Lacer is designed to close a medium sized orthosis or footwear size 42 EU and larger.

Comes pre assembled with a low profile dial and 8 guides. QuickFit Rivets sold separately, but are required to install the QuickFit Lacer kits. Lacers are sold individually. For a pair of shoes, two kits are required.

Technical data

Article number

514Z22=M-1



514Z22=M-10 QuickFit® Lacer

The medium size QuickFit Lacer is designed to close a medium sized orthosis or footwear size 42 EU and larger.

It comes with 10 Medium QuickFit Lacer kits and 100 QuickFit Rivets. For a pair of shoes, two kits are required.

Technical data

Article number

514Z22=M-10



514Z22=L-1 QuickFit® Lacer

Designed to close or secure above the ankle shoes or orthoses, such as a Leather AFO's, plastic rigid AFO's, or Crow boots. Can also be used on footwear such as hiking boots.

Comes pre assembled with a more powerful dial and 10 guides. Rivets sold separately.

Technical data

Article number

514Z22=L-1



21Z30=1 QuickFit® Releasing Lacer, LG

Adds opening functionality to our QuickFit Lacer by adding releasing guides to one side so that the system can fully open like a door for easy donning.

Comes with a more powerful reel to close larger devices. This kit is unassembled with 4 fixed guides and 3 releasing guides. Rivets sold separately.

Technical data

Article number

21Z30=1



21Z32=1 QuickFit® Tool

The QuickFit Tool is customized for use with QuickFit Rivets and is the recommended installation tool. The QuickFit Tool properly secures rivets and provides either an 1/8" and 3/16" rivet hole.

Technical data
Article number
21Z32=1



21Z31=1 QuickFit® Rivets 100-Pack

QuickFit Rivets come in packs of 100. You can supply your own rivets, or purchase our easy to use rivet pack to install QuickFit Lacers to any device or shoe.

Technical data
Article number
21Z31=1



21Z33=1 QuickFit® XS/S/M Dial and Lace Repl. Kit

Replacement BOA® dial and 1.5 m of Low-Mid Power Steel Lace for XS/SM/MD size QuickFit Lacers using BOA® dials.

Technical data
Article number
21Z33=1



21Z34=1 QuickFit® L Dial and Lace Repl. Kit

Replacement BOA® dial and 2 m of compatible Low-Mid Power Steel Lace for LG size QuickFit Lacers using BOA® dials.

Technical data
Article number
21Z34=1



21Z35=1 Low-Mid Power Steel Lace – 10 m Spool

Bulk replacement Low-Mid Power Steel Lace (10 m spool) compatible only with BOA® dials and QuickFit Lacer kits.

NOT COMPATIBLE with Click Reel.

Technical data

Article number

21Z35=1



514Z20=38-7 QuickFit strap

QuickFit Straps are designed to quickly install on a wide variety of cuffs, orthoses, and prosthetic devices to create an improved, easy to use, adjustable strap.

- Prefabricated, durable, and installs in minutes
- Don/doff via simple Velcro® closure
- 76mm (3") of adjustable travel
- Usable strap length: Max: 0.75m (30") / Min: 0.2m (9")
- Cut to size
- Width (38mm)



Information material

646D1441=EN_INT	514Z20/514Z21 Quick-Fit strap/QuickFit buckle produkt brief
-----------------	---

Technical data

Article number	Length	Width	Colour
514Z20=38-7	Approx. 73–63 cm	Approx. 38 mm	Black

Product example



BOA® is a registered trademark of Boa Technology, Inc.



514Z21=38-7 QuickFit buckle

The QuickFit Buckle can be added onto the QuickFit Strap in a few seconds. The magnetic buckle allows the user to easily connect the strap - even with limited dexterity or in hard to reach areas. Its plastic on plastic interface eliminates wear issues of Velcro® closure.

- Easily adds onto QuickFit Strap in seconds
- Adds 76mm (3") to strap usable length
- 5mm (.2") profile
- Width (38mm)

Technical data

Article number	Length	Width	Colour
514Z21=38-7	Approx. 18 cm	Approx. 38 mm	Black

Information material

646D1441=EN_INT	514Z20/514Z21 QuickFit strap/QuickFit buckle produkt brief
-----------------	--

Product example





21Z36=1 RevoFit® Re-Lacing Kit for BOA® dials

2 m of HD (High Durability) Lace for RevoFit Lamination Kits.
Compatible only with RevoFit Kits with BOA® dial.

Technical data

Article number

21Z36=1



21Z37=1 RevoFit® Dial and Lace Repl. Kit

Replacement BOA® dial and 2 m of HD (High Durability) Lace.
Compatible only with RevoFit Kits with BOA® dial.

Technical data

Article number

21Z37=1



21Z38=1 RevoLock® Lanyard Dial and Lace Kit

Replacement BOA® dial and 1.2 m of lace.

Compatible only with RevoLock Kits with BOA® dial. Not compatible with RevoLock V2 and V3 Kits

Technical data

Article number

21Z38=1



21Z39=1 RevoLock® Lanyard Re-Lacing Kit

Contains 1.2 m of lace.

Compatible only with RevoLock Kits with BOA® dial. Not compatible with RevoLock V2 and V3 Kits.

Technical data

Article number

21Z39=1



21Z40=1 **RevoLock® Lanyard – Threaded Insert**

Replacement RevoLock Lower-Extremity Threaded Insert which attaches the lace to the liner via a threaded connection.
Compatible only with RevoLock Lanyard and RevoLock Lanyard V2 Kits with BOA® dial.

Technical data
Article number
21Z40=1



21Z40=2 **RevoLock® Lanyard – Threaded Insert**

Replacement RevoLock Upper-Extremity Threaded Insert which attaches the lace to the liner via a threaded connection.
Compatible only with RevoLock Lanyard and RevoLock Lanyard V2 Kits with BOA® dial.

Technical data
Article number
21Z40=2

