

6Y112 *Skeo Sealing TT.*

The 6Y112 *Skeo Sealing TT* liner offers transtibial amputees a high level of comfort. The sealing ring reliably retains the vacuum in the socket. Its silky-smooth surface with no textile cover simplifies cleaning and putting on and taking off the prosthesis.

Designed for optimal outcomes

- A sealing ring reliably retains the vacuum in the socket
- Integrated sealing ring for the entire lifetime
- The inside of the liner is textured for a pleasant feeling on the skin
- Silky-smooth outer coating for easier donning and doffing
- Orientation aid on the inside of the liner for an optimal liner position on the residual limb
- An integrated matrix for reduced elongation and good tissue management
- The *Skeo Sealing TT* can be combined with the 21Y2 *QuickValve* or the *Harmony* system plus 21Y2 *QuickValve*.
- Protected against fresh, salt and chlorinated water



ottobock.



6Y112

Skeo Sealing TT.

Article number	Circumference	Ring height
6Y112=160X80	160 mm	80 mm
6Y112=180X80	180 mm	80 mm
6Y112=200X80	200 mm	80 mm
6Y112=210X80	210 mm	80 mm
6Y112=220X80	220 mm	80 mm
6Y112=235X80	235 mm	80 mm
6Y112=250X80	250 mm	80 mm
6Y112=265X80	265 mm	80 mm
6Y112=280X80	280 mm	80 mm
6Y112=300X80	300 mm	80 mm
6Y112=320X80	320 mm	80 mm
6Y112=340X90	340 mm	90 mm
6Y112=360X90	360 mm	90 mm
6Y112=380X90	380 mm	80 mm
6Y112=400X90	400 mm	90 mm



The ring



The material



The orientation aid