ottobock.



With BALLA

"

My prosthetist found a solution that allowed me to quickly return to my normal life.

László is an owner of a building company. After he lost his leg in an accident, for him it was essential to quickly return back to his normal life - with a prosthesis he could trust.

László was lucky and was fitted with a prosthesis that gave him a great deal of safety and independence right from the start. The adjustable socket in particular makes his daily routine easier. Just a few twists of the adjustment wheels and the socket fits. This also allows him to enjoy even longer walks with his dogs. László is relieved: "In the beginning, I could only walk for limited distances, but now I'm on the go from morning to evening again."

László's individual prosthetic solution.

Independent, safe, comfortable **The innovative socket solution**

Varos is a socket that benefits transfemoral amputees starting their rehabilitation phase as well as later during treatment with the final prosthesis. The socket circumference can be easily adjusted at any time using two adjustment wheels. This is important because the circumference of the residual limb changes, for example due to differences in temperature and also depending on how long it is worn during the day. Various activities are also easier to manage if the socket is somewhat snugger or looser.

Furthermore, the *Varos* socket and its *Varos* liner feature especially easy handling. The prosthesis can be easily put on and take off even while sitting; therefore the socket shells can be bent slightly.

Compared with a conventional socket, the *Varos* socket can be quickly custom fitted with only a few visits to a prosthetist.





Varos Liner

- Enables suspension with socket wall
- Improves pressure distribution

5A60 Varos socket

- Easy adaption of socket circumference by the user with BOA Closure System
- Length adaption due to cuttable shells & beams
- Flexible shells enable adaptivity
- Spheric adapter allows up to max. 15° Flexion (with additional adapters)



4R57 Rotation adapter

- Easy to put on and take off, and easy on your back
- Relaxed sitting

3C60 Kenevo

- Safety with every step
- Custom adaptation to the user's needs



1C11 Terion K2

- Safe heel strike
- Comfortable rollover and good energy return

For more information about Voros please visit.

#WeEmpowerPeople www.ottobock.com