ottobock.

Adapt. The mobility base.

The mobility base *Adapt* offers a wide range of adjustment options.

Thanks to the adjustable center of gravity, optimal driving behavior as well as an active and independent movement of the child can be encouraged.

Moreover, due to options such as drum brakes, anti-tippers, angle-adjustable arm supports or push bar, the mobility base *Adapt* provides easy handling and is an ideal companion in the everyday life.

Highlights

- For the child's independent and active mobility
- Seat tilt up to +50°
- Available in three sizes
- Variety of adjustment options





ottobock.

Adapt.The mobility base.



Available in three sizes

 The Adapt is available in three sizes and allows to adapt the seat width up to 10 cm without any additional parts.



Easy to use brakes

 The Adapt can be equipped with an drum brake or combi-drum brake system.



Functional seat tilt

 The Adapt has a functional seat tilt ranging from -10° up to +50° which allows different positioning.



Angle-adjustable push bar

• The push bar is angle-adjustable and provides thereby a comfortable handling.



Reliable anti-tippers

 The anti-tippers support safe driving and maneuvering of the mobility base.



Variety of adjustment options

 Due to the large variety of adjustment options, optimal driving characteristics can be achieved from active to passive mobility.



Practical luggage-bin

• Due to the practical luggage-bin, small items can be transported easily.



Angle-adjustable arm support

• The angle-adjustable arm support allows an individual adjustment to the user.



Technical Data

Feature	Adapt
Max. user weight	75 kg
Max. Payload	100 kg
Seat width range	Size 1: 30 – 40 cm Size 2: 40 – 50 cm Size 3: 50 – 60 cm
Base plate lenght	Size A: 32 cm Size B: 37 cm Size C: 42 cm
Tilt	-10 to +50°
Drive wheel sizes	20, 22, 24 and 26"
Crash tested (acc. ISO 7176-19)	✓
Weight incl. wheels	17,4 kg

Frame colors glossy (RAL-colors)



Frame colors matt (RAL-colors)

