

# 4R11 Quickchange. One turn – many options.





With the 4R11 *Quickchange* adapter, children and users with mobility grades 1-4 can remove their distal prosthetic components from the socket by themselves when needed, in just one step. This makes dressing and undressing easier and sitting more comfortable. In addition, a fitting with various feet or knee-foot combinations is possible, such as sports components.

## Product features

- Flexible fittings thanks to the four connection options – pyramid, pyramid receiver, thread and thread for modular system for children.
- The base element for all user groups is characterised by quick and easy handling. The locking mechanism can be operated with just one turn. The lock and arrow symbols allow the user to easily see in which direction the adapter should be opened and closed.
- With the 4X907=1, the pyramid receiver of the base element can be reduced so that the Quickchange adapter can also be used with the modular system for children.
- Protected against fresh, salt and chlorinated water.



## Technical data

				
Article# adapter	4R11=P-1-1	4R11=R-1-1	4R11=ST-1-1	4R11=ST-2-1
Article# base	4R11=R-3-2	4R11=R-3-2	4R11=R-3-2	4R11=R-3-2
Field of application	Adults/children: TT, KD, TF, HD (only adults)			
Mobility grade	1-4			
Max. body weight	150 kg			
Max. body weight for sports	125 kg			
Build height	44 mm	55 mm	41 mm	41 mm / 43 mm with 4X907=1
Weight	370 g	370 g	330 g	325 g
Connection	Pyramid/ pyramid receiver	Pyramid receiver/ pyramid receiver	Thread/ pyramid receiver	Thread for modular system for children/ pyramid receiver

## 4R11. Order form.

Complete and send order to: [order@ottobock.se](mailto:order@ottobock.se)







### Customer

Customer no.	
Company	
Street	
Postal code	City
O&P professional	

### Shipping address (if different from customer address)

Customer no.	
Company	
Street	
Postal code	City
Customer ref.	

### Single parts and accessories

Quantity	Article no.	Description	
	4R11=R-3-2	Base element with pyramid receiver	
	4R11=P-1-1	Connecting element with pyramid	
	4R11=R-1-1	Connecting element with pyramid receiver	
	4R11=ST-1-1	Connecting element with thread	
	4R11=ST-2-1	Connecting element with thread for modular system for children	
	4X907=1	Reduction pyramid receiver for modular system for children	

### Application examples

Easy removal of the prosthetic components below the socket for comfortable sitting/  
to change clothes more easily (1 set)



Change from everyday to sports prosthesis (1 set + 1 additional base element)



Change of two sockets (1 set + 1 additional connecting element)

