

4R11 Quickchange.

Amplify the possibilities.

With the 4R11 *Quickchange* adapter, children and users with mobility grades 1-4 can remove their distal prosthetic components from the socket by themselves when needed, in just one step. This makes dressing and undressing easier and sitting more comfortable. In addition, a fitting with various feet or knee-foot combinations is possible, such as sports components.

Product features

- Flexible fittings thanks to the four connection options – pyramid, pyramid receiver, thread and thread for modular system for children.
- The base element for all user groups is characterised by quick and easy handling. The locking mechanism can be operated with just one turn. The lock and arrow symbols allow the user to easily see in which direction the adapter should be opened and closed.
- With the 4X907=1, the pyramid receiver of the base element can be reduced so that the Quickchange adapter can also be used with the modular system for children.
- Protected against fresh, salt and chlorinated water.







Quick and easy handling

Four connection options

For adults and children

Technical data

				
Article number	4R11=P-1-3	4R11=R-1-3	4R11=ST-1-3	4R11=ST-2-3
Consists of	4R11=P-1-1 4R11=R-3-2	4R11=R-1-1 4R11=R-3-2	4R11=ST-1-1 4R11=R-3-2	4R11=ST-2-1 4R11=R-3-2
Field of application	Adults/children: TT, KD, TF, HD (only adults)			
Mobility grade	1-4			
Max. body weight	150 kg			
Max. body weight for sports	125 kg			
Build height	44 mm	55 mm	41 mm	41 mm / 43 mm with 4X907=1
Weight	370 g	370 g	330 g	325 g
Connection	Pyramid/ pyramid receiver	Pyramid receiver/ pyramid receiver	Thread/ pyramid receiver	Thread for modular system for children/ pyramid receiver

4R11. Order form.





Customer

Customer no.	
Company	
Street	
Postal code	City
O&P professional	







Shipping address (if different from customer address)

Customer no.	
Company	
Street	
Postal code	City
Customer ref.	

Set – consisting of 1 connecting element and 1 base element

Article no.	Image	Description	Quantity
4R11=P-1-3		Quickchange set with pyramid	
4R11=R-1-3		Quickchange set with pyramid receiver	
4R11=ST-1-3		Quickchange set with thread	
4R11=ST-2-3		Quickchange set with thread for modular system for children	

Single parts and accessories

Article no.	Image	Description	Quantity
4R11=R-3-2		Base element with pyramid receiver	
4R11=P-1-1		Connecting element with pyramid	
4R11=R-1-1		Connecting element with pyramid receiver	
4R11=ST-1-1		Connecting element with thread	
4R11=ST-2-1		Connecting element with thread for modular system for children	
4X907=1		Reduction pyramid receiver for modular system for children	

Application examples

Easy removal of the prosthetic components below the socket for comfortable sitting/to change clothes easier (1 set)



Order example

1x 4R11=P-1-3 Quickchange set with pyramid

Change from everyday to sports prosthesis and changing feet and knees
(1 set + 1 additional base element)



Order example

1x 4R11=P-1-3 Quickchange set with pyramid
1x base element 4R11=R-3-2

Change of two sockets
(1 set + 1 additional connecting element)



Order example

1x 4R11=P-1-3 Quickchange set with pyramid
1x connecting element with pyramid
4R11=P-1-1