

1C59 Taleo Adapt. Ready for everyday life.

The *Taleo Adapt* was designed for active users who navigate varied indoor and outdoor environments and place a high value on effortless walking wherever life takes them. Its hydraulic ankle unit allows the foot to adapt to uneven terrain and slopes for greater flexibility and comfort.

Designed for optimal outcomes

- Hydraulic ankle allows up to 2 degrees dorsiflexion and 10 degrees plantarflexion for a total range of 12 degrees.
- The dorsiflexion and plantarflexion adjustment valves allow resistance to be optimised to the individual user needs.
- Dorsiflexed position during swing phase for improved toe clearance.
- Three different angle indicator lines show the location of the hydraulic ankle joint (PF, neutral position, DF). These provide the technician with an orientation aid and enable a functionally correct alignment.
- Three different heel wedge options can be used to customise the impact stress at heel strike and the gait dynamics.
- Openings in the sole of the footshell prevent water from collecting in the prosthesis.
- Protected against fresh, salt and chlorinated water.

More than a foot. A foundation.



Mobility grade	3, 4
Max. body weight	130 kg
Side	left (L), right (R)
Sizes	22–30 cm
Weight without footshell*	622 g
Footshell colour	beige (4), light brown (15)
Weight with footshell*	899 g
Build height*	134 mm

* Reference size is 26 cm.



Hydraulic ankle
Adaptation to uneven surfaces
Improved socket comfort

1C59. Order form.

Customer

Customer no.	
Company	
Street	
Postal code	City
O&P professional	

Shipping address (if different from customer address)

Customer no.	
Company	
Street	
Postal code	City
Customer ref.	

Selection of the spring stiffness relative to

1. body weight and activity as well as

Body weight (kg)	Normal activity level	High activity level
Up to 51	1	2
52 to 58	2	3
59 to 67	3	4
68 to 77	4	5
78 to 88	5	6
89 to 100	6	7
101 to 115	7	8
116 to 130	8	8

2. foot size

Stiffness	Foot size (cm)							
	1	2	3	4	5	6	7	8
22						-	-	-
23							-	-
24							-	-
25								-
26								
27	-	-						
28	-	-						
29	-	-	-					
30	-	-	-					

Example for ordering 1C59=R-26-4-P/4N

Article no.	Side	Size	Stiffness	Colour	Shape
1C59	=	-	-	P /	N

Side	Size (cm)	Colour
R Right	22, 23, ..., 30	4 Beige
L Left		15 Light brown



The scope of delivery includes the 1C59 foot module, 2C8 footshell including connection cover, 2F50 heel wedge set (comprising three different degrees of hardness) and a black Spectra sock.

Recommended combinations*

TT	TF
21Y2=M10 QuickValve + 453A30/40 ProFlex Plus sleeve + 6Y95 Caleo 3D 21Y2=M10 QuickValve + 6Y112 Skeo Sealing TT	7T451=1 SiOCX TF Pro 6Y110 Skeo Sealing + 21Y21 Click Valve 21Y21 ClickValve + 6Y110 Skeo Sealing TF / 6Y111 Skeo Sealing 3D TF
	3C88-3=* / 3C98-3=* C-Leg 3B1-3=* Genium 3B5-3=* Genium X3 3B5-4=* Genium X4
4R11=* Quickchange	4R11=* Quickchange 4R57=WR Waterproof Rotation Adapter
	3F1=* + 99B120=* Functional cosmesis

* All components are sold separately and are available Ottobock products that are compatible with the 1C59 Taleo Adapt foot, which helps ensure optimal performance. Practitioners need to select components based upon individual patient criteria.