

C-Brace Interim KAFO.

Order form.


The **C-Brace Interim** fitting consists of industrially manufactured components that may be exclusively used on one single user for a maximum of 6 months. Longer use is not permitted. The **C-Brace** loaner joint must be ordered separately (see selection box below).

Contact person		Customer number		Date	
Customer					
Customer number					
Company					
Street					
Postcode		Town/City			
Shipping address (if different from customer address)					
Customer number					
Company					
Street					
Postcode		Town/City			
User data (optional)					
Patient ID					
Indication					
Weight		Side Left leg Right leg Both legs			

The following components are required for the fabrication of an interim orthosis:

C-Brace Interim KAFO

For thigh shell and AFO selection, please use the larger size in case of uncertainties.



Thigh shell

Side	Circumference A (size)	How to measure correctly	Article no.
left	30–45 cm (S)	Circumference measured 15 cm above knee pivot point (Nietert)	170K202=L-S
left	40–55 cm (M)		170K202=L-M
left	50–65 cm (L)		170K202=L-L
right	30–45 cm (S)		170K202=R-S
right	40–55 cm (M)		170K202=R-M
right	50–65 cm (L)		170K202=R-L

Knee section

Side	Variant (frontal knee angle)	Article no.	L	R
neutral	0°	Decision made by CPO during assembly	170K201=N-0	
neutral	3°		170K201=N-3	
neutral	6°		170K201=N-6	
neutral	9°		170K201=N-9	

AFO incl. ankle joint

Side	Ankle-bottom (B) / shoe size (C) (EU)	How to measure correctly	Article no.
left	70 mm / 36–41	Selection according to ankle-bottom measurement and shoe size	170K200=L-S
left	85 mm / 39–46		170K200=L-L
right	70 mm / 36–41		170K200=R-S
right	85 mm / 39–46		170K200=R-L

Accessories (reusable)

Article no.
KAFO Bag
642C3

C-Brace loaner joint

Side
17KO1=L-S
17KO1=R-S

Tools (reusable)

Article no.
C-Brace Interim Tool-Set
17OK203=N-1
C-Brace Interim Template-Set
17OK203=N-2