

Ottobock Kenevo Rehab Kit. Order form.

Complete the order form and send it to: order@ottobock.se

Customer

Company	
Customer Number	
Contact person	
Street	
Postal code	City
Email	
Order Number	

Weight		kg
Height		cm
Knee-floor measurement		cm
Date of amputation		
Side of amputation	☐ left ☐ right ☐ bilateral	

Amputation Level

- ☐ Knee Disarticulation
- ☐ Transfemoral
- ☐ Hip Disarticulation / Hemipelvectomy

Prosthetic attachment to residual limb

- ☐ Conventional socket
- ☐ Varos socket 5A60 (see separate order information)
- ☐ Osseointegration/bone anchored prosthesis (BAP)

Expected Mobility Grade

- ☐ 1 
- ☐ 2 
- ☐ 3 
- ☐ 4 
- ☐ not defined

Date of fitting

Knee joint (incl. AC adapter, inductive battery charger and AXON Tube Adapter)

QTY	Article No.	Description
	3C60=S	Kenevo
	3C60=ST-S	Kenevo for long residual limb
	2R17=___-S	AXON Tube Adapter

Foot

Article No.	Side	Size	–	Stiffness
☐ 1C11=			–	–S
☐ 1C30-1=			–	–S

The Kenevo Rehab Kit is provided for a maximum loan period of 6 months. The loaner fee will be charged after delivery of the Kenevo Rehab Kit according to the price list. If the Kenevo Rehab Kit is not returned in due time, a weekly loaner fee will be charged according to the price list. In addition, the general loaner conditions apply.

Shipping address (if different from customer address)

Company	
Street	
Postal code	City
Email	
Commission	

Stiffness chart Terion 1C11

Weight	Sizes			
	S (22–23 cm)	M (24–25 cm)	L (26–28 cm)	XL (29–30 cm)
up to 55 kg	1	1	2	–
56 – 75 kg	2	2	2	3
76 – 100 kg	3	3	3	3
101 – 125 kg	–	4	4	4
125 – 150 kg	–	5	5	5
151 – 175 kg	–	–	6	6

Stiffness chart Trias 1C30-1

Weight	Sizes					
	21–22 cm	23–24 cm	25–26 cm	27 cm	28 cm	29–30 cm
45 – 60 kg	1	1	–	–	–	–
61 – 80 kg	2	2	1	1	1	–
81 – 95 kg	–	3	2	2	2	1
96 – 110 kg	–	–	3	3	3	2
111 – 125 kg	–	–	–	4	4	3

Date

Place

Signature

Ottobock Kenevo Rehab Kit.

Order form.

Complete the order form and send it to: order@ottobock.se

Optional:

The Varos Socket System is not part of the rental Kenevo Rehab Kit but can be purchased additionally. Due to its volume adaptable features Varos is a perfect combination for the Kenevo Rehab Kit.

Varos Socket (Example 5A60=R-MC)

QTY	Article No.	Side (L , R)	–	Size (S, SC, M, MC, L)	Description
	5A60=		–		Varos Socket

Size	Varos S	Varos SC (conical)	Varos M	Varos MC (conical)	Varos L
Maximum body weight including additional weight	100 kg	100 kg	100 kg	100 kg	125 kg
Proximal residual limb circumference	410 – 510 mm	460 – 580 mm	460 – 580 mm	520 – 650 mm	520 – 650 mm
Distal residual limb circumference	320 – 390 mm	320 – 390 mm	350 – 440 mm	350 – 440 mm	400 – 490 mm
Residual limb length	200 – 320 mm	200 – 320 mm	200 – 320 mm	200 – 320 mm	200 – 320 mm

Varos Liner (Example 6Y200=M-2)

QTY	Article No.	Size	Description
	6Y200=		Varos Liner

Size	Proximal residual limb circumference	Distal residual limb circumference
S-1	420 – 490 mm	310 – 350 mm
SC-1 (conical)	480 – 550 mm	310 – 350 mm
S-2	450 – 520 mm	330 – 370 mm
SC-2 (conical)	510 – 580 mm	330 – 370 mm
M-1	480 – 550 mm	360 – 410 mm
MC-1 (conical)	530 – 610 mm	360 – 410 mm
M-2	510 – 580 mm	390 – 440 mm
MC-2 (conical)	560 – 650 mm	390 – 440 mm
L-1	530 – 610 mm	400 – 450 mm
L-2	560 – 650 mm	430 – 490 mm

 Please note, it is not possible to combine size L liners with sockets of any other size.

Date

Place

Signature