

# Triton 1C60

## Loaner Request Form

Email form to: [service\\_au@ottobock.com](mailto:service_au@ottobock.com)

Note: All loaners are subject to availability

ottobock.

Reference number:

Office use only

\* Required Fields

### Product Trial

Required Date\*

#### Trial Units

Ottobock Australia offers 1 week free patient trials, in addition to 1 week travel time to and from Sydney (2 weeks total).

A loaner fee of \$110 per week applies for each week the unit is outstanding after this time.

Patient Name / ID\*

Facility Name \*

Clinician's Name \*

Contact / Email \*

Delivery Address \*

If trial units are required for more than 1 week, please advise.

Number of weeks required

Purchase order number  
(if required)

### Example for ordering 1C60=R26-4-P/0-S

Article no.	Side	Size	- Stiffness	- P / 0-S
1C60=				- P / 0-S
1C60=				- P / 0-S

Side
Right R
Left L

Size [cm]
21, 22, ..... 30

### ① Selection Chart

Size	21	22	23	24	25	26	27	28	29	30
Body Weight KG										
up to 55	1	1	1	1	1	1	-	-	-	-
56 - 75	2	2	2	2	2	2	2	2	2	2
76 - 100	3	3	3	3	3	3	3	3	3	3
101 - 125	-	-	-	-	4	4	4	4	4**	4**
126 - 150	-	-	-	-	5	5	5*	5*	5*	5*

Slim footshell available  
(15 ± 5 mm heel height)

Both footshells available

Normal footshell available  
(10 ± 5 mm heel height)

\* When combining this configuration with the C-leg, please contact Ottobock Customer Service  
\*\* When combining this configuration with the C-leg or C-leg Compact, please contact Ottobock Customer Service

