

Triton Low Profile. Your will. Your way.

Triton Low Profile was designed for highly active individuals who navigate varied indoor and outdoor environments and place a high value on uncompromised response and control even when performing high-impact activities.

Designed for optimal outcomes

- An interconnected, linear spring design provides the level of midstance support needed for activities where responsiveness is essential
- For users with limited clearance
- Robust titanium adapter
- Customizable shock absorption with two different heel wedge options
- Suitable for a broad range of application, from everyday life to demanding professions and even leisure sports
- Slim footshell option
- Protected against fresh, salt, and chlorinated water

More than a foot. A foundation.



Technical data

Mobility grade	3, 4
Max. bodyweight	150 kg
Side	left (L), right (R)
Sizes	21–30 cm
Weight without footshell*	415 g
Footshell shape	Slim shape (S) with 15 ± 5 mm heel height (size 21–27 cm) Normal shape (N) with 10 ± 5 mm heel height. (size 24–30 cm)
Footshell colour	beige (4), light brown (15)
Weight with normal footshell*	635 g
Build height with normal footshell*	63 mm

* Reference size is 26 cm.



Strong and dynamic rollover
Energetic walking
Supports agile and high-impact movements

1C63

Order form.

Customer

Customer no.	
Company	
Street	
Postal code	City
O&P professional	

Shipping address (if different from customer address)

Customer no.	
Company	
Street	
Postal code	City
Customer ref.	

Selection of the spring stiffness relative to body weight and foot size

Foot size \ Body weight	Foot size										
	21 cm	22 cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm	
up to 55 kg	1	1	1	1	1	1	-	-	-	-	
56 – 75 kg	2	2	2	2	2	2	2	2	2	2	
76 – 100 kg	3	3	3	3	3	3	3	3	3	3	
101 – 125 kg	-	-	-	-	4	4	4	4	4*	4*	
126 – 150 kg	-	-	-	-	5	5	5*	5*	5*	5*	

- Slim footshell available (15 ± 5 mm heel height)
- Both footshells available
- Normal footshell available (10 ± 5 mm heel height)

* Do not combine this configuration with a 3C88-3/3C98-3 C-Leg 4.

Example for ordering 1C63=R27-3-P/4N

Article no.	Side	Size	Stiffness	Colour	Shape
1C63	=		-	P /	
	Side	Size (cm)	Colour	Shape	
	R Right	21, 22, ..., 30	4 Beige	S Slim 21-27 cm	
	L Left		15 Light brown	N Normal 24-30 cm	



The scope of delivery includes the 1C63 foot module, the 2C6 footshell including connection cover, the 2F60 heel wedge set (comprising two different degrees of hardness), and a black Spectra sock.

Recommended combinations**

TT	TF
 6Y95 Caleo 3D + 453A30/40 ProFlex Plus sealing sleeve + 21Y2=M10 QuickValve or 4R180 Harmony P4 + 6Y512=*-P Uneo Fresh + 453A30/40 ProFlex Plus sealing sleeve or 6Y112 Skeo Sealing TT + 21Y2=M10 QuickValve	 7T451=1 SiOCX TF Pro + 21Y21 ClickValve or 6Y110/6Y111 Skeo Sealing + 21Y2=M10 QuickValve
	 3B5-4 Genium X4 3C88-4/3C98-4 C-Leg 3R85 Dynion
 4R11 Quickchange	 4WR57=WR Rotation adapter 4R11 Quickchange
	 3F2=0 + 99B122* Functional cosmesis or 4P110=7/4P100 Protector long/short

** All components are sold separately and are available Ottobock products that are compatible with the 1C63 **Triton Low Profile** foot, which helps ensure optimal performance. Practitioners need to select components based upon individual patient criteria.