

Ottobock Kenevo Rehab Kit. Order form.

Complete the order form and send it to: order@ottobock.se

Customer

Company	
Customer Number	
Contact person	
Street	
Postal code	City
Email	
Order Number	

Shipping address (if different from customer address)

Company	
Street	
Postal code	City
Email	
Commission	

Weight kg

Height cm

Knee-floor measurement cm

Date of amputation

Side of amputation left right bilateral

Amputation Level

- Knee Disarticulation
- Transfemoral
- Hip Disarticulation / Hemipelvectomy

Prosthetic attachment to residual limb

- Conventional socket
- Varos socket 5A60 (see separate order information)
- Osseointegration/bone anchored prosthesis (BAP)

Expected Mobility Grade

- 1 
- 2 
- 3 
- 4 
- not defined

Date of fitting

Knee joint (incl. AC adapter, inductive battery charger and AXON Tube Adapter)

QTY	Article No.	Description
	3C60=3-S	Kenevo
	3C60=ST-3-S	Kenevo for long residual limb
	2R17=___-S	AXON Tube Adapter

Foot

Article No.	Side	Size	-	Stiffness
<input type="radio"/> 1C11=			-	-S
<input type="radio"/> 1C30-1=			-	-S

The Kenevo Rehab Kit is provided for a maximum loan period of 3 months. The loaner fee will be charged after delivery of the Kenevo Rehab Kit according to the price list. If the Kenevo Rehab Kit is not returned in due time, a weekly loaner fee will be charged according to the price list. In addition, the general loaner conditions apply.

Date Place Signature

Otto Bock Scandinavia AB · Postal: Box 4041, 169 04 Solna, Visiting: Barks väg 7, Solna · Sweden
T+ 46 11 280689 · info@ottobock.se · www.ottobock.se

Stiffness chart Terion 1C11

Weight	Sizes			
	S (22–23 cm)	M (24–25 cm)	L (26–28 cm)	XL (29–30 cm)
up to 55 kg	1	1	2	–
56 –75 kg	2	2	2	3
76 –100 kg	3	3	3	3
101 –125 kg	–	4	4	4
125 –150 kg	–	5	5	5
151 –175 kg	–	–	6	6

Stiffness chart Trias 1C30-1

Weight	Sizes					
	21–22 cm	23–24 cm	25–26 cm	27 cm	28 cm	29–30 cm
45 – 60 kg	1	1	–	–	–	–
61 – 80 kg	2	2	1	1	1	–
81 – 95 kg	–	3	2	2	2	1
96 – 110 kg	–	–	3	3	3	2
111 – 125 kg	–	–	–	4	4	3

Ottobock Kenevo Rehab Kit. Order form.

Complete the order form and send it to: order@ottobock.se

Optional:

The Varos Socket System is not part of the rental Kenevo Rehab Kit but can be purchased additionally. Due to its volume adaptable features Varos is a perfect combination for the Kenevo Rehab Kit.

Varos Socket (Example 5A60=R-MC)

QTY	Article No.	Side (L , R)	-	Size (S, SC, M, MC, L)	Description	
	5A60=		-		Varos Socket	
Size	Varos S		Varos SC (conical)	Varos M	Varos MC (conical)	Varos L
Maximum body weight including additional weight	100 kg		100 kg	100 kg	100 kg	125 kg
Proximal residual limb circumference	410 – 510 mm		460 – 580 mm	460 – 580 mm	520 – 650 mm	520 – 650 mm
Distal residual limb circumference	320 – 390 mm		320 – 390 mm	350 – 440 mm	350 – 440 mm	400 – 490 mm
Residual limb length	200 – 320 mm		200 – 320 mm	200 – 320 mm	200 – 320 mm	200 – 320 mm

Varos Liner (Example 6Y200=M-2)

QTY	Article No.	Size	Description
	6Y200=		Varos Liner
Size	Proximal residual limb circumference		Distal residual limb circumference
S-1	420 – 490 mm		310 – 350 mm
SC-1 (conical)	480 – 550 mm		310 – 350 mm
S-2	450 – 520 mm		330 – 370 mm
SC-2 (conical)	510 – 580 mm		330 – 370 mm
M-1	480 – 550 mm		360 – 410 mm
MC-1 (conical)	530 – 610 mm		360 – 410 mm
M-2	510 – 580 mm		390 – 440 mm
MC-2 (conical)	560 – 650 mm		390 – 440 mm
L-1	530 – 610 mm		400 – 450 mm
L-2	560 – 650 mm		430 – 490 mm

i Please note, it is not possible to combine size L liners with sockets of any other size.

Responsible physiotherapist:

Date

Place

Signature