

Taleo Harmony.

Ready for everyday life.

Taleo Harmony was designed for active individuals who navigate varied indoor and outdoor environments and place a high value on effortless walking and the ability to go wherever life takes them.

Designed for optimal outcomes

- Integrated Harmony P3 pump improves control over the prosthesis, minimises stump problems and skin irritation and compensates for fluctuations in volume
- Effectively absorbs torsion to relieve the residual limb and boost comfort in everyday life in combination with vertical shock absorption
- Dual springs and a long carbon base enable a smooth rollover and efficient energy return at varying speeds
- Thanks to three different heel wedges, both the impact stress at heel strike and the gait dynamics can be customised to individual needs

More than a foot.
A foundation.



Technical data

Mobility grade	2, 3, 4
Max. body weight	150 kg
Side	left (L), right (R)
Sizes	22 – 30 cm
Weight without footshell*	751 g
Footshell shape	Slim shape (S) with 15 +/- 5 mm heel height (size 22 – 25 cm) Normal shape (N) with 10 +/- 5 mm heel height (size 24 – 30 cm)
Footshell colour	beige (4), light brown (15)
Weight with normal footshell*	980 g
Build height with normal footshell*	185 mm
Other features	Vacuum; +/- 10° torsion movement; up to 15 mm vertical shock absorption

* Reference size is 26 cm.



Firm hold and control throughout the day

Relief of residual limb and enhanced comfort

Smooth rollover

Energy efficient walking

Terrain conformance on uneven ground

1C52

Order form.

Customer

Customer no.	
Company	
Street	
Postal code	City
O&P professional	

Shipping address (if different from customer address)

Customer no.	
Company	
Street	
Postal code	City
Customer ref.	

Selection of the spring stiffness relative to

1. body weight and activity as well as

Body weight (kg)	Normal activity level	High activity level
up to 51	1	2
52 – 58	2	3
59 – 67	3	4
68 – 77	4	5
78 – 88	5	6
89 – 100	6	7
101 – 115	7	8
116 – 130	8	9
131 – 150	9	–

2. foot size

Slim footshell available (15 ± 5 mm heel height)
Both footshell available
Normal footshell available (10 ± 5 mm heel height)

Stiffness	Foot size (cm)								
	1	2	3	4	5	6	7	8	9
22						–	–	–	–
23							–	–	–
24							–	–	–
25								–	–
26								*	–
27	–	–						*	*
28	–	–						*	*
29	–	–	–					*	*
30	–	–	–					*	*

* Do not combine this configuration with a 3C88-3/3C98-3 C-Leg 4.

Example for ordering 1C52=R26-4-P/4N

Article no.	Side	Size	Stiffness	Colour	Shape
1C52	=		-	- P /	
	Side	Size (cm)	Stiffness	Colour	Shape
	R Right	22, 23, ..., 30	1, 2, ..., 9	4 Beige	S Slim 22-25 cm
	L Left			15 Light brown	N Normal 24-30 cm



The scope of delivery includes the foot module 1C52, the 2C15 footshell including the connection cap, the 2F50 heel wedge set (comprising three different degrees of hardness), a 2Z360 functional ring set for replacement, a 2R117=0 socket connector, a pre-compression kit and a black Spectra sock. Please order the 4Y383 cosmetic exhaust flange for a cosmetic cover separately.

Recommended combinations**

TT



453A30/40=* ProFlex Plus knee sleeve + 6Y512=*-P Uneo Fresh
or
453A30/40=* ProFlex Plus knee sleeve + 6Y400 Uneo Unique



Aluminium line adapters
4R11 Quickchange

** All components are sold separately and are available Ottobock products that are compatible with the 1C52 **Taleo Harmony** foot, which helps ensure optimal performance. Practitioners need to select components based upon individual patient criteria.