

## Release Note | CG000547 | Rev A

# Chromium X Series (X/iX) Firmware Release Note

For use with:

**Chromium X** & Accessory Kit, 12-Month Warranty, PN-1000331 Chromium iX Accessory Kit, PN-1000323 | Chromium X Upgrade Kit, PN-1000327

**Chromium X** & Accessory Kit, 24-Month Warranty, PN-1000332 Chromium iX Accessory Kit, PN-1000323 | Chromium X Upgrade Kit, PN-1000327

**Chromium iX** & Accessory Kit, 12-Month Warranty, PN-1000328 Chromium iX Accessory Kit, PN-1000323

**Chromium iX** & Accessory Kit, 24-Month Warranty, PN-1000329 Chromium iX Accessory Kit, PN-1000323

## Notices

#### **Document Number**

CG000547 | Rev A

#### **Legal Notices**

© 2022 10x Genomics, Inc. (10x Genomics). All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of 10x Genomics, is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. 10x Genomics provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. A non-exhaustive list of 10x Genomics' marks, many of which are registered in the United States and other countries can be viewed at: www.10xgenomics.com/trademarks. 10x Genomics may refer to the products or services offered by other companies by their brand name or company name solely for clarity, and does not claim any rights in those third-party marks or names. 10x Genomics products may be covered by one or more of the patents as indicated at: www.10xgenomics.com/patents. The use of products described herein is subject to 10x Genomics Terms and Conditions of Sale, available at www.10xgenomics.com/legal-notices, or such other terms that have been agreed to in writing between 10x Genomics and user. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.Notice to Customer: The listed products and their use are the subject of United States Patent Nos. 6,159,736, and 6,294,385, European Patent No. 1115856 and related patents and patent applications licensed from the Wisconsin Alumni Research Foundation.

#### **Instrument & Licensed Software Updates Warranties**

Updates to existing Instruments and Licensed Software may be required to enable customers to use new or existing products. In the event of an Instrument failure resulting from an update, such failed Instrument will be replaced or repaired in accordance with the 10x Limited Warranty, Assurance Plan or service agreement, only if such Instrument is covered by any of the foregoing at the time of such failure. Instruments not covered under a current 10x Limited Warranty, Assurance Plan or service agreement will not be replaced or repaired.

#### Support

Email: **support@10xgenomics.com** 10x Genomics 6230 Stoneridge Mall Road Pleasanton, CA 94588 USA

### Chromium X Series Firmware Release Note

	Standards
For Use With	PN-1000331 Chromium X & Accessory Kit, 12-Month Warranty PN-1000332 Chromium X & Accessory Kit, 24-Month Warranty PN-1000328 Chromium iX & Accessory Kit, 12-Month Warranty PN-1000329 Chromium iX & Accessory Kit, 24-Month Warranty
Firmware Version	1.1.0
Release Date	March 2022
Compatibility	Chromium X Series
Features	<ul> <li>New assay compatibility- Fixed RNA Profiling</li> <li>Enhanced error detection and troubleshooting capabilities</li> <li>Ethernet and Wi-Fi MAC addresses are accessible</li> <li>The instrument's serial number is set as the hostname enabling all Chromium X Series instruments to have a unique hostname</li> </ul>
Bug Fixes	<ul><li>Fixed Wi-Fi icon display in the connection status</li><li>Disabled bypassing critical components installation during firmware update</li></ul>