




CryoPro® Handling devices




Carrying bag

Ideal for transporting your CryoPro® when working remotely. The bag has a practical strap for safety in car.




Protective gloves and eyewear

We provide protective wear for both hands and eyes when filling your CryoPro® bottle.



Filters and gaskets

Replaceable filters for withdrawal tube, and gaskets which ensure the tightness of your CryoPro® bottle.




Roller base

Universal type for easy move of dewars with liquid nitrogen. Fits different dewar sizes.



Withdrawal tube

Universal type for pouring the liquid nitrogen from your dewar to your CryoPro® bottle.



Storage containers

We offer a wide selection of metal dewars in different sizes and static holding time.

CryoPro®

A versatile low-cost treatment
of common skin lesions

Cortex

Niels Jernes Vej 6B
9220 Aalborg
Denmark
+45 9857 4100
cortex@cortex.dk
www.cortex.dk



CryoPro® Treatment accessories

The CryoPro® device is delivered in a box with some of our most common treatment and handling accessories as a standard

Spray Tips



Spray Tip "A"
1.00 mm

*In the CryoPro® box



Spray Tip "B"
0.75 mm

*In the CryoPro® box



Spray Tip "C"
0.55 mm (pre-mounted)


*In the CryoPro® box




Spray Tip "D"
0.45 mm

*In the CryoPro® box


Cones



Lexan Disc
with 4 cones




Rubber Cones
set with 6 sizes



Acne Soft
Spray Tip



Straight Spray
extension, 0.55 mm



Bent Spray
extension, 0.55 mm

*In the CryoPro® box

Attachments



Spray Tip
cleaning adapter

*In the CryoPro® box




Luer Lock
adapter



Open Flow
restrictor

Contact Probes



Sharp pointed
Contact probe



1 mm
Contact probe



2 mm
Contact probe



3 mm
Contact probe



4 mm
Contact probe



5 mm
Contact probe



6 mm
Contact probe




8 mm
Contact probe




10 mm
Contact probe



15 mm
Contact probe



20 mm
Contact probe



30 mm
Contact probe

As cool as it gets –

CryoPro® delivers safety, reliability,
and precision at -196°C



Cortex

Niels Jernes Vej 6B
9220 Aalborg
Denmark
+45 9857 4100
cortex@cortex.dk
www.cortex.dk