




CryoPro® Handling devices




Carrying bag

Ideal for transporting your CryoPro® when working remotely. The bag has a practical strap for safety in car.




Protective gloves and eyewear

We provide protective wear for both hands and eyes when filling your CryoPro® bottle.



Filters and gaskets

Replaceable filters for withdrawal tube, and gaskets which ensure the tightness of your CryoPro® bottle.




Roller base

Universal type for easy move of dewars with liquid nitrogen. Fits different dewar sizes.



Withdrawal tube

Universal type for pouring the liquid nitrogen from your dewar to your CryoPro® bottle.



Storage containers

We offer a wide selection of metal dewars in different sizes and static holding time.

CryoPro®

A versatile low-cost treatment
of common skin lesions

Cortex

Niels Jernes Vej 6B
9220 Aalborg
Denmark
+45 9857 4100
cortex@cortex.dk
www.cortex.dk




CryoPro® Treatment accessories

The CryoPro® device is delivered in a box with some of our most common treatment and handling accessories as a standard

We offer a complete range of cryosurgery equipment including treatment devices and accessories for cryosurgery.


We can provide closed contact probes, open spray tips, extensions, adapters and spray cones.

Spray Tips




Spray Tip "A" 1.00 mm

*In the CryoPro® box




Spray Tip "B" 0.75 mm

*In the CryoPro® box




Spray Tip "C" 0.55 mm (pre-mounted)

*In the CryoPro® box




Spray Tip "D" 0.45 mm


*In the CryoPro® box




Conical Spray Tip "A" 1.00 mm



Conical Spray Tip "B" 0.75 mm



Conical Spray Tip "C" 0.55 mm




Conical Spray Tip "D" 0.45 mm

Attachments




Replaceable Nitrogen Filter (pre-mounted)

*In the CryoPro® box



Spray Tip cleaning adapter




Luer Lock adapter




Open Flow restrictor

Cones



Lexan Disc with 4 cones



Rubber Cones set with 6 sizes



Straight Spray extension "A" 1.00 mm, 25 mm short



Straight Spray extension "B" 0.75 mm, 25 mm short



Straight Spray extension "C" 0.55 mm, 25 mm short



Straight Spray extension "C" 0.55 mm, 75 mm long




Bent Spray extension "B" 0.75 mm, 75 mm long



Bent Spray extension "C" 0.55 mm, 75 mm long




Acne Soft Spray Tip



Stand for Spray Tips

Contact Probes



Sharp pointed Contact probe



1 mm Contact probe



2 mm Contact probe



3 mm Contact probe



4 mm Contact probe



5 mm Contact probe



6 mm Contact probe



8 mm Contact probe



10 mm Contact probe



15 mm Contact probe



20 mm Contact probe



30 mm Contact probe

As cool as it gets –

CryoPro® delivers safety, reliability, and precision at -196°C

Cortex
Niels Jernes Vej 6B
9220 Aalborg
Denmark
+45 9857 4100
cortex@cortex.dk
www.cortex.dk