

Explore our
**DermaLab
Aesthetic**

True skin analysis for clinics

Dedicated to unparalleled
client satisfaction?

Prepared to upscale and
increase sales
conversions?

Revolutionize
your skincare business
with DermaLab Aesthetic!

The True Skin Analysis

The DermaLab Aesthetic is a powerful skin analysis tool that, within a few minutes, provides a comprehensive overview of the client's skin, covering up to five highly relevant parameters.

The results are communicated through user-friendly software, which also offers printouts/e-mails for the clients allowing them to remember and reference your services. This scientifically grounded tool ensures precise and valid data – the true skin analysis.

Attract clients – increase sales

Invite clients to your clinic by offering them the true skin analysis and providing thorough insight into their skin. Attract even hesitant and uncertain clients in search of reassurance by empowering them with data.


When your clients understand the true needs of the skin, the aesthetic skin treatment decisions become easy. Use the DermaLab Aesthetic to assist your clients in taking the next step and increase your sales conversions.

Optimize client satisfaction and loyalty

The skin analysis serves as a crucial starting point for treatment planning and facilitates optimal expectation management, ultimately leading to the best possible outcomes and client experiences.

Employ the DermaLab Aesthetic during follow-up consultations with clients to assess and demonstrate results and to plan further actions.

Skin analysis powerhouse
Understanding runs deep
for those who analyze

 **Cortex**
Niels Jernes Vej 6B
9220 Aalborg
Denmark
+45 9857 4100
cortex@cortex.dk
www.cortex.dk

DermaLab Aesthetic

First edition of our new tool for true skin analysis in clinics



Extensive skin experience

Behind this indispensable tool is the Danish company Cortex Technology with more than 30 years of experience in the production of precision skin tools and a trusted partner in research. All products are meticulously designed, manufactured, and assembled by hand in Denmark, with global distribution partnerships in place.

Essential parameters

The DermaLab Aesthetic equips you with five essential parameters:




The Ultrasound Probe gives you an instant assessment of the collagen level



The pH Probe measures skin surface pH



The Elasticity Probe measurement of the skin's ability to bounce back when it is stretched



The Hydration Probe measures the water binding capacity of the outermost layer of the skin



The Sebum Module provides a quick assessment of skin oiliness

Made by Cortex

Don't settle for a quick fix
Make the right choice easy

Cortex

Niels Jernes Vej 6B
9220 Aalborg
Denmark
+45 9857 4100
cortex@cortex.dk
www.cortex.dk