

| Description | Size | Part Number |
|----------------------|---------------|-------------|
| Switchcut Reamer Kit | 4.5 x 6.0 mm | 110027674 |
| | 4.5 x 6.5 mm | 110027675 |
| | 4.5 x 7.0 mm | 110027676 |
| | 4.5 x 7.5 mm | 110027677 |
| | 4.5 x 8.0 mm | 110027678 |
| | 4.5 x 8.5 mm | 110027679 |
| | 4.5 x 9.0 mm | 110027680 |
| | 4.5 x 9.5 mm | 110027681 |
| | 4.5 x 10 mm | 110027682 |
| | 6.0 x 11.0 mm | 110027684 |
| 6.0 x 12.0 mm | 110027686 | |

| Description | Part Number |
|-----------------------------------|-------------|
| SwitchCut Universal Guide Body | 110026899 |
| ACL Fem Guide Arm - Point | 110026880 |
| ACL Fem Guide Arm - U-Shape Right | 110026900 |
| ACL Fem Guide Arm - U-Shape Left | 110026901 |
| PCL Tibial Guide Arm | 110026881 |
| Meniscal Root Tibial Guide Arm | 110026882 |
| ACL Tibial Guide Arm - Hole | 110026905 |
| Tibial Guide Arm Point to Point | 110026903 |
| SwitchCut Guide Bullet 4.5 mm ID | 110026898 |
| SwitchCut Guide Bullet 6.0 mm ID | 110026902 |
| Centering Sleeve 2.4 mm | 110026904 |
| SwitchCut Caddy | 110026911 |

Introducing the Complete Switchcut™ Line Up



All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for health care professionals. Distribution to any other recipient is prohibited.

For indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the package insert or contact your local representative; visit www.zimmerbiomet.com for additional product information.

Not for distribution in France.

Check for country product clearances and reference product specific instructions for use.

©2020 Zimmer Biomet

Authorized Representative
Biomet Sports Medicine
56 East Bell Drive
P.O. Box 587
Warsaw, Indiana 46581
USA

Legal Manufacturer
T.A.G. Medical Products
Corporation LTD
Kibbutz Gaaton 2513000
Israel
Tel: +972-4-9858400
Fax: +972-4-9858404

www.zimmerbiomet.com

 **ZIMMER BIOMET**
Your progress. Our promise.®

2609.1 GLBL-en REV1220

 **ZIMMER BIOMET**

SwitchCut™ Reaming System

The SwitchCut Reaming System is multifunctional, easy to use, and intuitive. Used to create bone sockets and tunnels for ACL, PCL, and Meniscal Root procedures, the SwitchCut reamers automatically transition from 4.5 or 6.0mm tunnel diameter to the appropriate socket diameter by simply running the reamer in reverse. Never lose the bone tunnel because the SwitchCut reamers are cannulated allowing for simple relay suture passage in one step. Pick from six SwitchCut Guide Arms to mate with the Universal Guide Body and Bullet to accommodate many surgical techniques.

System consists of left or right specific guide, bullet, and reamer with Nitinol passer

- 6–12 mm reamer sizes, including half sizes
- Quick and easy deployment of reamer blade by simply switching power on the drill from forward to reverse
- Blade will deploy once it contacts bone in a reverse setting
- Cannulated reamers allow for passing of the Nitinol loop through the tunnel to assist in keeping tunnel placement



Femoral Guide Arm - Right

**Tibial Guide
Point to Point**

PCL Tibial Arm Guide

ACL Femoral Crescent Guide

Meniscal Root Tibial Guide Arm

ACL Fem Guide Arm - Point

