



# Persona<sup>®</sup>

REVISION KNEE SYSTEM



# THE POWER TO **PERSONALIZE**<sup>™</sup>

Persona Revision offers anatomic implant geometry to provide patients an individualized fit while achieving the desired bone coverage, fixation<sup>1</sup> and stability. The comprehensive instrumentation platform<sup>2</sup> empowers surgeons with the flexibility to tailor their surgical approach intraoperatively, while the streamlined delivery system's modularity is designed to help increase OR efficiencies and workflow. Combined with the clinically successful Trabecular Metal<sup>™</sup> Material<sup>3-5</sup> and Vivacit-E<sup>®</sup> Vitamin-E Highly Crosslinked Polyethylene, Persona Revision truly offers personalization for revision knee arthroplasty.





# PERSONALIZED IMPLANTS ACHIEVING AN INDIVIDUALIZED FIT.

## **Anatomic Implants**

The anatomic design of the Persona Primary Tibia was incorporated into the Persona Revision Tibial Implant to provide a unified system approach to address primary and revision knee arthroplasty.

## **Achieving Zonal Fixation**

One study shows that in order to achieve robust implant fixation, most successful knee revision procedures require a combination approach to fixation with multi-zone capabilities.<sup>6</sup> Persona Revision offers numerous anatomic tibial and femoral cones, augments and an array of stem options to address fixation needs.

## **Restoring Stability**

One of the most common causes of implant failure is joint instability.<sup>7-9</sup> In an effort to address this, Persona Revision's multiple femoral components in varying sizes, along with a vast continuum of bearing constraints provide surgeons the tools to help restore a balanced knee.

# PRECISE INSTRUMENTATION WITH A PERSONALIZED APPROACH.

## **Comprehensive Instrument Platform**

Persona Instrumentation was designed to be versatile in its capabilities and philosophies, precise in its measurements<sup>10-13</sup>, comprehensive, and comfortable with repetitive use.

## **Personalized Surgical Techniques**

The surgeon defined workflow allows surgeons to individualize the surgical experience based on the patient's bone and joint condition.

## **Streamlined Delivery System**

Persona Revision's efficient kitting strategy with modular brackets allows for adaptable tray configurations based on surgeon and hospital preference.





# PROVEN TECHNOLOGY


BUILT ON CLINICAL PERFORMANCE.

## **Trabecular Metal™ Technology**

This highly porous material is designed to replicate the structure, function and physiological properties of cancellous bone<sup>14</sup> with over 20 years of clinical use.<sup>15</sup>

## **Vivacit-E® Vitamin-E HXPE**

Vivacit-E Polyethylene is actively stabilized with Vitamin-E to help protect against oxidation and maintain wear resistance and strength throughout the life of the implant. Exceptional oxidative stability,<sup>16, 17</sup> ultra-low wear<sup>18, 19</sup> and enhanced strength.<sup>20,21</sup>



EXPANDING THE  
**Personalized**  
KNEE FAMILY

Persona Revision Knee



PERSONALIZED IMPLANTS.  
**PRECISE INSTRUMENTATION.**  
PROVEN TECHNOLOGY.



Persona Primary Knee



Persona Partial Knee

1. Zimmer ZTR\_BI\_0245-18.
2. Persona Revision Knee Surgical Technique. 1832-GLBL-en.
3. Zhang, Y., et al. Interfacial Frictional Behavior: Cancellous Bone, Cortical Bone, and a Novel Porous Tantalum Biomaterial. *Journal of Musculoskeletal Research*. 3(4):, 245-251, 1999.
4. Bobyn, J.D., et al. Characteristic of Bone In-growth and Interface Mechanics of a New Porous Tantalum Biomaterial. *Journal of Bone and Joint Surgery (British)*. 81-B(5):907, 1999.
5. Shirazi-Adl, A., et al. Experimental Determination of Friction Characteristics at the Trabecular Bone / Porous-coated Metal Interface in Cementless Implants. *The Journal of Biomedical Research*. 27: 167-175, 1993.
6. Morgan-Jones, J., et al. Zonal Fixation in Revision Total Knee Arthroplasty. *Bone and Joint Journal*. 97-B: 147-9, 2015.
7. Dalury, D., et al. Why Are Total Knee Arthroplasties Being Revised? *The Journal of Arthroplasty*. 28(8): 120-121, 2013.
8. Sharkey, P.F., et al. Why Are Total Knee Arthroplasties Failing Today-Has Anything Changed After 10 Years? *The Journal of Arthroplasty*. 29(7): 1774-1778, 2014.
9. Chang, M.J., et al. Diagnosis, Causes and Treatments of Instability Following Total Knee Arthroplasty. *Knee Surgery and Related Research*. 26(2): 61-67, 2014.
10. Internal Data on File at Zimmer Biomet. Zimmer ZRR\_WA\_2409\_11. 2011.
11. Internal Data on File at Zimmer Biomet. Zimmer ZRR\_WA\_2403\_11. 2011.
12. Internal Data on File at Zimmer Biomet. Zimmer ZRR\_WA\_2551\_12. 2012.
13. Internal Data on File at Zimmer Biomet. Zimmer ZTR\_BI\_0133\_18. 2018.
14. Bobyn, J.D., et al. Characterization of New Porous Tantalum Biomaterial for Reconstructive Orthopaedics. Scientific Exhibition. 66th Annual Meeting of the American Academy of Orthopaedic Surgeons. Anaheim, CA. 1999.
15. Levine, B. et al. Experimental and Clinical Performance of Porous Metal Tantalum in Orthopedic Surgery. *Biomaterials*. 27: 4671-81, 2006.
16. Zimmer ZRR\_WA\_2409\_11.
17. Zimmer ZRR\_WA\_2403\_11.
18. Zimmer ZRR\_WA\_2537\_12.
19. Zimmer ZTR\_BI\_0249\_18.
20. Zimmer ZRR\_WA\_2551\_12.
21. Zimmer ZTR\_BI\_0133\_18.

All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet.

This material is intended for healthcare professionals, Zimmer Biomet employees, the Zimmer Biomet sales force, and authorized representatives. Distribution to any other recipient is prohibited. For indications, contraindications, warnings, precautions, potential adverse effects and patient counselling information, see the package insert or contact your local representative; visit [www.zimmerbiomet.com](http://www.zimmerbiomet.com) for additional product information.

Results are not necessarily typical, indicative, or representative of all recipient patients. Results will vary due to health, weight, activity and other variables. Not all patients are candidates for this product and/or procedure. Only a medical professional can determine the treatment appropriate for your specific condition. Appropriate post-operative activities and pain will differ from patient to patient. Talk to your surgeon about whether joint replacement is right for you and the risks of the procedure, including the risks of infection, implant wear, loosening, breakage or failure, any of which can necessitate additional surgery or treatment. For additional information or to find a surgeon near you, visit [www.zimmerbiomet.com](http://www.zimmerbiomet.com) [or call 1-800-HIP-KNEE].

Not for distribution in France.

Check for country product clearances and reference product specific instructions for use.

©2019, 2022 Zimmer Biomet



2730.2-GLBL-en-Issue Date 2022-01-12



**Legal Manufacturer**

Zimmer, Inc.  
1800 W. Center Street  
Warsaw, Indiana 46580  
USA



CE mark on a surgical technique is not valid unless there is a CE mark on the product label.



**Authorized Representative**

Zimmer GmbH  
Sulzerallee 8  
8404 Winterthur  
Switzerland

[www.zimmerbiomet.com](http://www.zimmerbiomet.com)