

El Passo™ Cannula



Don't go with the flow

Flexible Silicone Cannula
designed for improved
visualization and maneuverability
within the repair site.¹

- **Double Dam** design helps manage fluid and prevent excessive leakage.
- The 5mm **adjustable** spacer creates different length options, allowing for one cannula to be used across different patient sizes.
- Low profile **distal flange** sits flush against tissue for improved visualization and reduced risk of unintentional pullout.
- **Flexible design** facilitates passing of implants and curved instruments for ease of maneuvering within the repair site.



Part Number	Size	Color
110048000	8 mm x 30 mm	Gray
110048001	8 mm x 40 mm	Gray
110048002	8 mm x 50 mm	Gray
110048004	10 mm x 50 mm	Gray

1. Improved visualization and maneuverability over rigid body cannula.

Confidential information intended for the Zimmer Biomet sales force or authorized representative. Distribution to any other recipient is prohibited. All content herein is protected by copyright, trademarks and other intellectual property rights, as applicable, owned by or licensed to Zimmer Biomet or its affiliates unless otherwise indicated, and must not be redistributed, duplicated or disclosed, in whole or in part, without the express written consent of Zimmer Biomet. For indications, contraindications, warnings, precautions, potential adverse effects and patient counseling information, see the Instructions for use or contact your local representative; visit www.zimmerbiomet.com for additional product information. ©2024 Zimmer Biomet

 **ZIMMER BIOMET**
Moving You Forward.™

4481.1-GLBL-en Issue Date 2024-07