

Product highlights

- Practical multifunctional table increasing patient comfort
- Can be installed on either side of seat
- Space-saving solution (no need for any additional mobile bed-side table)
- Easy and fast adjusting by hand-screw
- Can be tilted and folded aside for easy patient entry and leaving the chair
- Plastic tray is easy to remove and clean



Comfort-T2 multifunctional table
with tiltable tray



Technical data

Load capacity	10 kg
Size of tray	40 x 56 cm
Material	powder coated steel & plastic
Tilt angle of tray	0° / 15°

Sitting / Entry position

Adjustability of the integrated table and tilt-up armrest allow the patient to easily enter or leave the chair

