



AED Monthly Checklist

15 Feb 2021

Complete

Score	13 / 16 (81.25%)	Flagged items	1	Actions	2
Location of AED	2114 Central St 6th floor, Kansas City, MO 64108, USA (39.0867543, -94.5879528)				
AED Serial No.	12345678910				
Checked on	15 Feb 2021 07:15 PST				
Checked by	SafetyCulture Staff				

Flagged items & Actions

1 flagged, 2 actions

Flagged items

1 flagged, 0 actions

AED Monthly Checklist

Are the electrode pads in the AED container?

No



Photo 2

Other actions

2 actions

AED Monthly Checklist / Check PPE/Ready Kit

Scissors



To do | Priority: Medium | Due: 19 Feb 2021 08:00 PST | Created by: SafetyCulture Staff

Buy scissors and place in AED ready kit.

AED Monthly Checklist

Test AED container tamper alarm



In progress | Priority: High | Due: 15 Feb 2021 08:00 PST | Created by: SafetyCulture Staff

Test AED container tamper alarm.

AED Monthly Checklist

1 flagged, 2 actions, 12 / 15 (80%)

Is the AED secure in its container?

Yes



Photo 1

Check AED is on standby mode - green flashing light



If AED is not on standby mode, verify AED operation

Verify AED Operation

4 / 4 (100%)

Open AED lid / Turn on AED



Wait for status indicator to change from RED to GREEN after five seconds



Close AED lid / Turn off AED



Wait for status indicator to change from RED to GREEN after five seconds



Are the electrode pads in the AED container?

No



Photo 2

Check expiration date of electrode pads



Electrode Pads Expiration Date:

22 Feb 2021 00:00 PST

Is there an extra set of electrode pads?

Yes

Check expiration date of AED battery



The AED battery is expired if 4 years have passed since its installation.

AED Battery Installation Date:

14 Feb 2020 00:00 PST

Check PPE/Ready Kit

1 action, 3 / 4 (75%)

Mask



Scissors



To do | Priority: Medium | Due: 19 Feb 2021 08:00 PST | Created by: SafetyCulture Staff

Buy scissors and place in AED ready kit.

Razor



Gloves



Test AED container tamper alarm



In progress | Priority: High | Due: 15 Feb 2021 08:00 PST | Created by: SafetyCulture Staff

Test AED container tamper alarm.

IMPORTANT Maintenance

1 / 1 (100%)

If electrode pads were expired, did you replace them?

N/A

Electrode pads not expired.

Input the date the electrode pads were replaced:

16 Dec 2020 00:00 PST

If AED battery was expired, did you replace it?

N/A

AED battery not expired.

Input the date the AED battery was replaced:

14 Feb 2020 00:00 PST

If the AED container tamper alarm was not working, did you change its battery?

Yes



Photo 3

Signature

A handwritten signature in black ink, appearing to be 'Aed'.

AED Master
15 Feb 2021 07:41 PST

Media summary



Photo 1



Photo 2



Photo 3