



QUICK REFERENCE LAYERING GUIDE

	Dentin Layer	Enamel Layer	Optional Effect	Result
Bleach or Junior	B0.5 (A0.5 or B1)	Enamel White 	Enamel Trans at incisal edge 	
Young	B1 	Enamel Blush 	Highlights with Opaque White 	
Adult	A2 (A1 or A3)	Enamel Neutral 	—	
Adult Senior	A3 	Enamel Yellow 	—	
Posterior Single Shade	—	Enamel Neutral 	—	
Posterior Double Shade	A2 	Enamel Neutral 	—	