

where
STRONG
bonds
begin



opaletch™
phosphoric acid gel

Achieve the ideal bonding surface for complete resin penetration

35% PHOSPHORIC ACID

Opal Etch 35% phosphoric acid etchant is based on the award-winning, top-rated Ultra-Etch[™] etchant chemistry. The unique physical and chemical properties of Opal Etch etchant create the optimal viscosity to ensure it won't migrate or prematurely dry. Its syringe delivery with the Opal[™] Blue tip facilitates precise placement and complete control, and its self-limiting formulation creates the ideal surface for adhesive penetration.

Strong bonds begin with Opal Etch etchant.

- Optimal viscosity prevents migration or premature drying
- Proprietary etch solution rinses off enamel quickly and easily
- Brilliant blue color assures safe, accurate placement and complete removal
- Contains a surfactant for complete, easy rinsing
- Unique self-limiting properties prevent over-etching



500-012 - Opal Etch Syringe Kit
4 x 1.5 g Opal Etch syringes
20 x Opal Blue tips



Brilliant blue color makes Opal Etch etchant easy to see, while its optimal viscosity ensures no migration.



Deliver Opal Etch etchant with the Opal Blue tip for precise placement.



500-013 - Opal Etch Syringe 20pk
500-091 - Opal Etch Syringe 50pk
1.2 ml syringe



500-090 - Opal Etch IndiSpense Syringe 1pk
30 ml syringe

710 - Inspiral Brush Tips 20pk



124 - 1.2 ml Empty Delivery Syringe 20pk
1.2 ml syringe

500-019 - Opal Blue Tips 20pk

