

Complete protection
you can see



OpalDam® Green

Light-cured Resin Barrier

OpalDam Green resin barrier offers effective coverage of oral tissues, making in-office whitening easier than ever before.

ULTRADENT
PRODUCTS, INC.



800.552.5512 ULTRADENT.COM

OpalDam® Green
Light-cured Resin Barrier

OpalDam® Green Resin Barrier

- Vibrant green color is easily seen for complete coverage
- Light-reflective properties focus heat away from teeth and gingiva for patient comfort
- Viscous solution is easy to apply on tissue, resulting in fast preparation times

OpalDam Green methacrylate-based resin barrier is used to isolate tissue during in-office tooth whitening procedures. The green color is easy to see against gingiva and teeth, resulting in a complete barrier. OpalDam Green resin barrier is easy to apply and remove and has been specifically designed to easily break from undercuts during removal.



Viscous solution is easy to apply.



Light-reflective properties focus heat away from teeth and gingiva for patient comfort.



Complete barrier is ideal for tooth whitening.



OpalDam Green resin cleanly breaks from undercuts.

1824 OpalDam Green Kit

*4 x 1.2 ml syringes
10 x Micro 20 ga tips
10 x Black Mini tips*



OpalDam Green 1.2 ml Syringes

1825 4pk 1826 20pk



Refrigerate if not used on a daily basis.

OpalDam® Resin Barrier

- Passive adhesive (sealing)
- Methacrylate-based resin

OpalDam resin barrier is used to isolate tissue adjacent to teeth being whitened. For single-tooth whitening, OpalDam resin barrier may be used to protect adjacent teeth.



Viscous solution is easy to apply.



OpalDam resin cleanly breaks from undercuts.

324 OpalDam Kit

*4 x 1.2 ml syringes
10 x Micro 20 ga tips
10 x Black Mini tips*



OpalDam 1.2 ml Syringes

325 4pk 326 20pk



Refrigerate if not used on a daily basis.