



The Gemini 810 + 980 diode laser is dentistry's first dual-wavelength soft tissue diode laser featuring 20 watts of peak super-pulsed power and a stunning, fully transparent electroluminescent display. No matter the procedure, the innovative Gemini laser makes it faster, smoother, and more efficient.

DUAL WAVELENGTHS.
ONE OF A KIND.


GEMINI™
810 + 980 DIODE LASER



- 20 watts of peak super-pulsed power for faster, smoother cutting
- Dual wavelength technology combines the optimal melanin absorption of the 810 nm wavelength and the optimal water absorption of the 980 nm wavelength
- Tip illumination provides better visibility at surgical site
- Sleek, innovative design features a stunning transparent electroluminescent display
- Simple user interface and 19 preset procedures enhance ease of use
- Wireless foot pedal and battery operation allow for convenient movement from operator to operator
- Autoclavable handpiece for simple sterilization between procedures
- Designed and assembled in the USA from U.S. and imported components



Diode Laser Peak Power Comparison¹



8980 - Gemini Laser Kit

1 x Each power supply, foot pedal, and handpiece
3 x Safety glasses sets
10 x 5 mm tips

8998 - Gemini PBM Adapter Kit

1 x Photobiomodulation (PBM) adapter
2 x Spacers
1 x Handpiece holder clip
1 x Cleaning cloth

8981 - Gemini Power Supply 1pk

8983 - Gemini 5 mm Pre-Initiated Tip 25pk 8984 - Gemini 7 mm Uninitiated Tip 25pk

8999 - Gemini PBM Spacer Tip 5pk

8985 - Gemini Safety Glasses 1pk

8986 - Handpiece Shell 1pk

8982 - Gemini Foot Pedal 1pk

Note: If you would like to order a Gemini Foot Pedal, please contact Equipment Repair at (801) 553-4574.

*SIROLaser Advanced™ is a trademark of Sirona Dental; EPIC™ 10 is a trademark of Biolase, Inc.; Precise® is a registered trademark of CAO Group, Inc.; Picasso™ is a trademark of AMD Laser. 1. Data published by manufacturer 2. Peak power in dual wavelength mode.

