

DUAL WAVELENGTHS. ONE OF A KIND.



The Gemini 810 + 980 diode laser is dentistry's first dual-wavelength soft tissue diode laser featuring 20 watts of peak super-pulsed power and a stunning, fully transparent electroluminescent display. No matter the procedure, the innovative Gemini laser makes it faster, smoother, and more efficient.

GEMINI™

810 + 980 DIODE LASER



ULTRADENT
PRODUCTS, INC.


800.552.5512 ULTRADENT.COM



- 20 watts of peak super-pulsed power for faster, smoother cutting
- Dual wavelength technology combines the optimal melanin absorption of the 810nm wavelength and the optimal water absorption of the 980nm wavelength
- Tip illumination provides better visibility at surgical site
- Sleek, innovative design features a stunning transparent electroluminescent display
- Simple user interface and 19 preset procedures enhance ease of use
- Wireless foot pedal and battery operation allow for convenient movement from operator to operator
- Autoclavable handpiece for simple sterilization between procedures
- Designed and assembled in the USA from U.S. and imported components



Diode Laser Peak Power Comparison¹



8980 Gemini Laser Kit
1 x Power supply
1 x Foot pedal
1 x Handpiece
10 x 5 mm tips
3 x Safety glasses sets

8981 Gemini Power Supply

8982 Gemini Foot Pedal

8983 Gemini 5 mm Tip - 25pk

8984 Gemini 7 mm Tip - 25pk

8985 Gemini Safety Glasses

8986 Handpiece Shell

*SIROLaser Advanced™ is a trademark of Sirona Dental; EPIC™ 10 is a trademark of Biolase, Inc.; Precise® is a registered trademark of CAO Group, Inc.; Picasso™ is a trademark of AMD Laser.

1. Data published by manufacturer
2. Peak power in dual wavelength mode.