

VALO™ X

BROADBAND LED CURING LIGHT



NEW FINISHES
SAME IMPRESSIVE POWER

NEW FINISHES

- Now you can express your individual style in the office with a vibrant, timeless color

HIGH-GRADE AEROSPACE ALUMINUM

- Unibody design is exceptionally durable and allows for excellent thermal management

INCREASED LENS SIZE

- 12.5 mm lens covers any tooth while maintaining accessibility and patient comfort

ACCELEROMETER FUNCTION

- Allows you to quickly and easily move through curing and diagnostic aid modes

SIMPLIFIED INTERFACE

- Curing and diagnostic aid modes are indicated, operated, and activated with top and bottom buttons or by Accelerometer Function

12 LED CHIPSET

- Provides high-intensity, broadband light for excellent beam uniformity, curing depth, and beam collimation

MULTI-CONFIGURATION

- Can be used in a corded or cordless configuration

5-YEAR WARRANTY

- Includes a 5-year Manufacturer Warranty

CURING MODES: Standard Power Mode, Xtra Power Mode | **DIAGNOSTIC AID LIGHT MODES:** White Light Diagnostic Aid Mode, Black Light Diagnostic Aid Mode

CYCLING BETWEEN MODES



CURING MODES:
Move the VALO X light forward in a drum tap motion to access and cycle through curing modes.



DIAGNOSTIC AID LIGHT MODES:
Move the VALO X light in a drum tap motion to the side to access and cycle through diagnostic aid light modes.

LENS ACCESSORIES

TWO CURING LENSES



ProxiCure™ Ball Lens

5889



PointCure™ Lens

5890

THREE DIAGNOSTIC AID LENSES



TransLume™ Lens

5891



Diffuser Lens

5892



Interproximal Lens

5893

VALO X Colors Kit –

- 1 x VALO X LED curing light
- 1 x Diffuser lens
- 2 x Rechargeable batteries
- 1 x Battery charger
- 1 x Blue-light-blocking shield
- 1 x Sample pack of barrier sleeves
- 1 x Power supply (for battery charger or cord adapter*)

- 6496 Red
- 6497 Pink
- 6498 Gray
- 6499 Green
- 6519 Blue
- 6538 Purple
- 6539 Black

5973 - VALO X Kit –



- 1 x VALO X LED curing light
- 5 x Accessory lenses
- 2 x Rechargeable batteries
- 1 x Battery charger
- 1 x Power supply (for battery charger or cord adapter)
- 1 x Cord adapter
- 1 x Handpiece bracket holder
- 1 x Blue-light-blocking shield
- 1 x Sample pack of barrier sleeves



SCAN
TO VIEW
PRICING

Accessories

- 4951 VALO X Battery Charging Unit
- 4952 VALO X Power Supply
- 5189 VALO X Assembled Cord

- 5437 VALO X Batteries 2pk
- 4665 VALO X Barrier Sleeves 100pk
- 5898 VALO X Barrier Sleeves 600pk
- 5839 Curing Light Ease-In-Shields Kit

*Cord adapter not included.

© 2025 Ultradent Products, Inc. All rights reserved.
1011267-PKAR03 082525

800.552.5512 | ULTRADENT.COM