


Better tools for better results

At Ultradent, we appreciate your precise craftsmanship and know that your hard work enhances patient comfort and helps to improve oral health. Our quality line of lab products is designed to make your job easier from start to finish.



Whitening Tools

Dental Lab Products and Tools

Ultradent™ LC Block-Out Resin

Ultradent LC Block-Out Resin is adept at filling in voids and blocking out undercuts. With its high contrast, optimal viscosity, and syringe delivery, it is easy to place accurately and precisely. It is also ideal for creating reservoir space in whitening trays. Studies show that reservoirs can improve whitening results by keeping more gel in contact with teeth.

Provides reservoir space for whitening trays and is useful for other laboratory procedures such as model and die repairs.




240 - Ultradent LC Block-Out Resin 4pk Kit
4 x 1.2 ml syringes
20 x Black Mini tips



241 - Ultradent LC Block-Out Resin Syringe 4pk
1.2 ml syringes

242 - Ultradent LC Block-Out Resin 20pk Kit
20 x 1.2 ml syringes
20 x Black Mini tips



Ultradent LC Block-Out Resin is a hard, strong, no-mix material for blocking out undercuts on dies and filling in voids.

Sof-Tray™ Classic Sheets

Select the 0.035" for most whitening trays, and the 0.060" or the 0.080" for whitening patients who are bruxers.

226 - Sof-Tray Regular Sheets 0.035" 5" x 5" 25pk
0.9 mm - 127 x 127 mm
0.035" thickness



227 - Sof-Tray Medium Sheets 0.060" 5" x 5" 20pk
0.5 mm - 127 x 127 mm
0.060" thickness

284 - Sof-Tray Heavy Sheets 0.080" 5" x 5" 20pk
2.0 mm - 127 x 127 mm
0.080" thickness

Ultradent™ Utility Vinyl Cutters

Ideal for gross trimming of tray.



604 - Ultradent Utility Vinyl Cutters 1pk

Ultradent™ Ultra-Trim Scalloping Scissors

- For precision trimming of border around interdental papilla
- Spring loaded to minimize finger fatigue
- Grip tray material easily
- Made of durable stainless steel



605 - Ultradent Ultra-Trim Scalloping Scissors 1pk

