

NEW!



GEMINITMNOVA

810 + 980 DIODE LASER



TOMORROW'S DENTISTRY IN A CORDLESS LASER

The Gemini NovaTM laser is the only cordless soft tissue diode laser that combines 20 watts of peak super-pulsed power for faster and gentler cutting, dual wavelengths for clinical versatility, and numerous benefits of photobiomodulation, all in a compact, ultra portable solution that fits seamlessly into your practice.^{1,2}

Why choose the Gemini Nova laser?



The cordless wand and compact base deliver true cord-free operation, eliminating fiber hassles and improving the overall user experience

20 WATTS

The 20 watts of peak, super-pulsed power provide faster cutting, less heat, and ultra clean incisions in soft tissue¹



Three wavelength modes



Take full advantage of laser benefits with the three included photobiomodulation adapters: 3 mm, 7 mm, and 25 mm



Adjust power wirelessly with a simple twist of the wand battery



Preset procedures are divided into intuitive categories for efficient use



HD capacitive touch display offers easy access to preset procedure settings, providing the best user experience yet



Custom-designed wand and base user-replaceable battery packs eliminate downtime and service interruptions



Wi-Fi connectivity allows for over-the-air updates and dedicated tech support



Mobile app and dashboard monitor usage statistics, including ROI and procedure data



Uses the same disposable fiber tips as the previous Gemini™ lasers



Two-year warranty included



6351 - Gemini Nova Kit

- 1 x Gemini Nova laser base
- 1 x Wireless laser wand
- 2 x Wand battery packs
- 1 x Foot pedal
- 2 x Safety glasses
- 10 x 5 mm Disposable fiber tips
- 1 x Power supply
- 1 x Base battery pack
- 3 x PBM adapters (25 mm, 7 mm, 3 mm)
- 5 x PBM light adapters
- 1 x Barrier sleeve starter pack



8983 - Gemini 5 mm Pre-Initiated Tip 25pk

8984 - Gemini 7 mm Uninitiated Tip 25pk



5837 - Dual Wave Ease-In-Shields Kit 1pk

OD 5+ @ 800-900nm OD 6+ @ 900-980nm
OD 7+ @ 980-1064nm OD 6+ @ >1064-1080nm



8985 - Gemini Safety Glasses 1pk

Note: The items above are shared between Gemini, Gemini EVO, and Gemini Nova lasers.



6362 - Gemini Nova 25 mm PBM Spacer Tips 5pk

6361 - Gemini Nova PBM Light Adapters 10pk



6355 - Gemini Nova Wand Battery 1pk



6364 - Gemini Nova 3 mm PBM Adapter Kit

1 x 3 mm PBM adapter, 2 x PBM light adapters

6366 - Gemini Nova 7 mm PBM Adapter Kit

1 x 7 mm PBM adapter, 2 x PBM light adapters

6363- Gemini Nova 25 mm PBM Adapter Kit

1 x 25 mm PBM adapter, 2 x PBM light adapters
2 x 25 mm PBM spacers



6358 - Gemini Nova Base Battery Pack 1pk



6356 - Gemini Nova Power Supply - U.S. 1pk

6360 - Gemini Nova Fiber Tip Barrier Sleeve 100pk
6914 - Gemini Nova 3 mm PBM Tip Barrier Sleeve 50pk
6915 - Gemini Nova 7 mm PBM Tip Barrier Sleeve 50pk