


# THE ORIGINAL **KNITTED CORD**


Easy packing. Excellent absorption. Exceptional retraction.

# Ultrapak™ Knitted Cord

When **tissue management is key**, reach for Ultrapak knitted displacement cord. With easy **packing**, excellent **absorption**, and exceptional **retention**, the proprietary knitted design of Ultrapak cord has been the preferred choice for years. As an important part of Ultradent's proven tissue management system, Ultrapak cord displaces tissue with predictability and minimum effort.



Knitted Ultrapak cord is composed of thousands of tiny, interlocking loops so it compresses and expands easier than other cords. 100% cotton fibers provide high absorption<sup>1</sup> of hemostatic agents and sulcular fluids. Use cords that appear slightly larger than you are used to with braided cords.



Ultrapak cord packs effortlessly into the sulcus and then exerts a gentle, continuous outward force as knitted loops seek to open for optimal tissue displacement.



Ultrapak CleanCut design features a blade in the cap for efficient cutting. A special dispensing orifice prevents cord from falling into the bottle.

Ask about our Mix and Match program!

137 - Ultrapak Cord #000 Regular 1pk



132 - Ultrapak Cord #1 Regular 1pk



199 - Ultrapak Cord #1 Epinephrine 1pk

136 - Ultrapak Cord #00 Regular 1pk



133 - Ultrapak Cord #2 Regular 1pk



197 - Ultrapak Cord #00 Epinephrine 1pk

200 - Ultrapak Cord #2 Epinephrine 1pk

131 - Ultrapak Cord #0 Regular 1pk



134 - Ultrapak Cord #3 Regular 1pk



198 - Ultrapak Cord #0 Epinephrine 1pk

**130 - Ultrapak Kit**  
1 x Each #00, 0, 1, and 2 cord  
1 x Ultrapak organizer



**Recommended Product: Fischer's Ultrapak Packers**  
Specially designed packers ease the packing of Ultrapak knitted cord.  
Available with thin serrated or non-serrated edges.



SCAN  
TO VIEW  
PRICING

1. Data on file.

© 2023 Ultradent Products, Inc. All rights reserved.  
87577AR13 102522

800.552.5512 | ULTRADENT.COM

