


Complete protection you can see



OpalDam™ Green

Light-Cured Resin Barrier

OpalDam Green resin barrier offers effective coverage of oral tissues, making in-office whitening easier than ever before.



OpalDam™ Green
Light-Cured Resin Barrier

OpalDam™ Green Resin Barrier

- Vibrant green color is easily seen for complete coverage
- Light-reflective properties focus heat away from teeth and gingiva for patient comfort
- Viscous solution is easy to apply on tissue, resulting in fast preparation times

OpalDam Green methacrylate-based resin barrier is used to isolate tissue during in-office tooth whitening procedures. The green color is easy to see against gingiva and teeth, resulting in a complete barrier. OpalDam Green resin barrier is easy to apply and remove and has been specifically designed to easily break from undercuts during removal.

Refrigerate if not used on a daily basis.



Viscous solution is easy to apply.



Light-reflective properties focus heat away from teeth and gingiva for patient comfort.



Complete barrier is ideal for tooth whitening.



OpalDam Green resin cleanly breaks from undercuts.



1824-U - OpalDam Green Syringe Kit
 4 x 1.2 ml OpalDam Green syringes
 10 x Black Mini tips
 10 x Micro 20 ga tips



1825-U - OpalDam Green Syringe 4pk
1826-U - OpalDam Green Syringe 20pk
 1.2 ml syringes

OpalDam™ Resin Barrier

- Passive adhesive (sealing)
- Methacrylate-based resin

OpalDam resin barrier is used to isolate tissue adjacent to teeth being whitened. For single-tooth whitening, OpalDam resin barrier may be used to protect adjacent teeth.

Refrigerate if not used on a daily basis.



Viscous solution is easy to apply.



OpalDam resin cleanly breaks from undercuts.



324-U - OpalDam Syringe Kit
 4 x 1.2 ml OpalDam syringes
 10 x Black Mini tips
 10 x Micro 20 ga tips



325-U - OpalDam Syringe 4pk
326-U - OpalDam Syringe 20pk
 1.2 ml syringes

