



### Peeso Reamer Bur

- Prepares canal orifice for straight line access
- Can be used for initial removal of gutta-percha in retreatment
- Non-cutting tip for increased safety



Refill and Assorted 6pk

No.	Tip Size	19 mm
1	0.50 mm	3729A
2	0.70 mm	3730A
3	0.90 mm	3731A
4	1.10 mm	3732A
5	1.30 mm	3733A
6	1.50 mm	3734A
Asst	1-6	3728A

### Gates Glidden Bur

- Allows for deeper penetration in straight canals
- Can also be used in straight part of curved canals
- Non-cutting tip for increased safety



Refill and Assorted 6pk

No.	Tip Size	19 mm
1	0.70 mm	3736A
2	0.90 mm	3737A
3	1.10 mm	3738A
4	1.30 mm	3739A
5	1.50 mm	3740A
6	1.70 mm	3741A
Asst	1-6	3735A

## Genius™

reciprocation & rotary motor

- Get the safety of reciprocation (90°/30°) and the efficiency of rotary (360°)
- Complete any procedure with just two files
- Multiple finishing options for final enlargement
- Designed for treatment and retreatment

ALSO  
AVAILABLE  
FROM  
ULTRADENT

## Genius™

reciprocation & rotary motor



Use Genius Files with the Genius Endodontic motor and contra-angle

30	.08	3743A	
		18 mm	
25	.04	3746A	3747A
		21 mm	25 mm
30	.04 JJ	3748A	3749A
		21 mm	25 mm
35	.04 J	3750A	3751A
		21 mm	25 mm
40	.04 I	3754A	3755A
		21 mm	25 mm
50	.04 I	3756A	3757A
		21 mm	25 mm



# ENDO-EZE™

Ultradent Endodontics

## INTRODUCING ENDO-EZE FILES

Endo-Eze stainless steel and nickel-titanium files are created using a state-of-art alloy treatment, giving them the flexibility to negotiate curvatures and the resistance to advance to the apical third. The ergonomic handles improve clinical performance, giving better tactile feedback and reducing fatigue.

Simply, our files are designed to be the best!

### Stainless Steel K-Files

- Balanced flexibility and resistance
- 60° angle non-cutting tip for effortless file penetration
- Spiral design with four cutting edges (square cross section)



Refill 6pk

Color	Tip Size	21 mm	25 mm	31 mm
	08	3619A	3620A	3621A
	10	3622A	3623A	3624A
	15	3625A	3632A	3633A
	20	3634A	3635A	3638A
	25	3639A	3640A	3641A
	30	3642A	3643A	3644A
	35	3645A	3646A	3647A
	40	3648A	3649A	3653A
	45	3654A	3655A	3656A
	50	3657A	3658A	3659A
	55	3660A	3661A	3663A
	60	3664A	3665A	3666A
	70	3667A	3668A	3669A
	80	3670A	3671A	3672A

Assorted 6pk

Color	Tip Size	21 mm	25 mm	31 mm
Asst	15-40	3673A	3675A	3677A
Asst	45-80	3674A	3676A	3678A

Refill 6pk

Color	Tip Size	21 mm	25 mm
	15	4256A	4258A
	20	4257A	4259A
	25	3611A	3615A
	30	3612A	3616A
	35	3613A	3617A
	40	3614A	3618A

Assorted 6pk

Color	Tip Size	21 mm	25 mm
Asst	15-40	3609A	3610A

### NiTi K-Files

- Nickel-titanium thermal treated alloy
- Provides sharpness and cutting ability with added flexibility
- Ideal for enlargement of severely curved canals



### Hedstrom (H-Files)

- Spiral designed flutes increase in depth and length from tip to handle
- Enhanced efficiency while maintaining safety
- Specially designed for retreatment



Refill 6pk

Color	Tip Size	21 mm	25 mm	31 mm
	15	3686A	3692A	3698A
	20	3687A	3693A	3699A
	25	3688A	3694A	3700A
	30	3689A	3695A	3711A
	35	3690A	3696A	3712A
	40	3691A	3697A	3713A

Assorted 6pk

Color	Tip Size	21 mm	25 mm	31 mm
Asst	15-40	3683A	3684A	3685A

### Flexfiles

- Stainless steel alloy with high flexibility
- Square cross-section design with more cutting edges than K-files
- Non-cutting tip for increased safety



Refill 6pk


Color	Tip Size	21 mm	25 mm
	15	3716A	3722A
	20	3717A	3723A
	25	3718A	3724A
	30	3719A	3725A
	35	3720A	3726A
	40	3721A	3727A

Assorted 6pk

Color	Tip Size	21 mm	25 mm
Asst	15-40	3714A	3715A

### Finger Spreaders

- Safely opens space for lateral compaction of gutta-percha
- Tapered design
- ISO standard sizing for gutta-percha compatibility



Assorted 6pk

Color	Tip Size	25 mm
Asst	20-30	3742A