

1 > 2

Say goodbye to multiple band cements.

Opal Band Cement's unique triple cure formula is versatile enough to work for both regular and fixed appliance banding—*so one cement is all you'll ever need.*

opal[™] band cement



opal™ band cement



opal™ band cement is everything you're looking for in an orthodontic band cement. Its versatile chemistry is strong enough for fixed appliance forces, yet easier to remove than a resin modified glass ionomer, and its convenient delivery system makes application a breeze!

Versatile Chemistry

Opal Band Cement eliminates the need to have one cement for regular banding and one cement for fixed appliance banding. Its unique triple cure formula creates optimal bond strengths for both! Opal Band cement also contains glass ionomer fillers that react with moisture to strengthen the cure and release fluoride.

Convenient Delivery

Say goodbye to mixing powders and liquids. Opal Band Cement comes in a convenient dual-barrel syringe that mixes the chemistry for you. It creates the ideal viscosity for easy placement and removal, and makes application simple and mess-free.

Available in 4pk syringe kits.



500-000 - Opal Band Cement 4pk Syringe Kit
4 x 8.5 g Opal Band syringes
50 x Mixing tips

500-086 - Mixing Tip Refill Kit
50 x Mixing tips



Opal Band Cement's bright blue chemistry is easy to see whether you're applying or removing it.



SCAN
TO VIEW
PRICING



© 2023 Ultradent Products, Inc. All rights reserved.
62277AR06 103122

800.552.5512 | ULTRADENT.COM

